# ElAgen CEA kit

REF LI4022F1

$\Sigma$  96 tests


**Adaltis Italia S.p.A.**  
Via Cristoni, 12  
40033 Casalecchio di Reno – (BO) Italy  
Tel. +39-051-6136511 – Fax + 39-051-575280  
[www.adaltis.com](http://www.adaltis.com)

FOR IN VITRO DIAGNOSTIC USE ONLY 

Store at 2...8 °C

en

## SYMBOLS USED ON LABELS

<b>MTPLATE</b>	Microplate
<b>CONJ HRP</b>	HRP Conjugate
<b>CAL0</b>	Calibrator 0
<b>CAL1</b>	Calibrator 1
<b>CAL2</b>	Calibrator 2
<b>CAL3</b>	Calibrator 3
<b>CAL4</b>	Calibrator 4
<b>CAL5</b>	Calibrator 5
<b>SUBS TMB</b>	TMB substrate
<b>SOLN STOP</b>	Stop solution
<b>WASH BUF 50X</b>	Washing buffer concentrate
<b>LOT</b>	Lot number
<b>REF</b>	Catalogue Code
	Expiry date (Use by...)
 2°C 8°C	Temperature limitation (store at 2...8°C)
$\Sigma$	Number of tests
	Keep away from sunlight
	Manufactured by
	Attention, See Instructions For Use
<b>IVD</b>	In vitro diagnostic medical device (In vitro diagnostic use)
	Biological risks (see SECTION 7)- Warnings and Precautions).

### 1.0 INTENDED USE

ElAgen CEA kit is a direct solid phase enzyme immunoassay for quantitative determination of Carcinoembryonic Antigen (CEA) in human serum.

### 2.0 SUMMARY AND EXPLANATION OF THE TEST

Carcinoembryonic antigen (CEA) is a glycoprotein with a molecular weight of 180 kDa. CEA is the first of the so-called carcinoembryonic proteins that was discovered in 1965 by Gold and Freeman (1).

CEA is the most widely used marker for gastrointestinal cancer. Although CEA is primarily associated with colorectal cancers, other malignancies that can cause elevated levels of CEA include breast, lung, stomach, pancreas, ovary and other organs. Benign conditions that cause significantly higher than normal levels include inflammation of lung and gastrointestinal (GI) tract and benign liver cancer (2,3).

Heavy Smokers, as a group, have higher than normal baseline concentration of CEA.


### 3.0 PRINCIPLE OF THE TEST

#### Immunoenzymometric assay:

In this method, calibrators, patient specimens and/or controls (containing the native CEA antigen) are first added to streptavidin coated wells. Biotinylated monoclonal and enzyme labeled antibodies are then added and the reactants mixed: these antibodies have high affinity and specificity and are directed against distinct and different epitopes of CEA.

Reaction between the various CEA antibodies and native CEA occurs in the microwells without competition or steric hindrance, forming a soluble sandwich complex

The interaction is illustrated by the following equation:

**BtnAb(m)** = Biotinylated Monoclonal Antibody (Excess Quantity)

**AgCEA** = Native Antigen (Variable Quantity)


**EnzAb** = Enzyme labeled Antibody (Excess Quantity)

**EnzAb - AgCEA - BtnAb(m)** = Antigen-Antibodies Sandwich Complex

**k<sub>a</sub>** = Rate Constant of Association

**k<sub>-a</sub>** = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

**StreptavidinC.W.** = Streptavidin immobilized on well

**Immobilized complex** = sandwich complex bound to the well

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is quantitated by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

### 4.0 CONTENT OF THE KIT

#### 4.1 MICROPLATE (code 66PIA)

1 microplate of 12 strips x 8 wells coated with Streptavidin. and packaged in aluminium bag with a drying agent.  
Store at 2-8°C.

#### 4.2 CALIBRATORS

6 vials of 1ml each containing serum calibrators for CEA in human serum. Preservative has been added.

Ready to Use

Store at 2-8°C.

Calibrator	symbol	code	Concentration ng/ml	Concentration (mIU/ml)
Calibrator 0	<b>CAL0</b>	66CAL0	0 ng/ml	0 mIU/ml
Calibrator 1	<b>CAL1</b>	66CAL1	5 ng /ml	50 mIU/ml
Calibrator 2	<b>CAL2</b>	66CAL2	10 ng /ml	100 mIU/ml
Calibrator 3	<b>CAL3</b>	66CAL3	25 ng /ml	250 mIU/ml
Calibrator 4	<b>CAL4</b>	66CAL4	50 ng /ml	500 mIU/ml
Calibrator 5	<b>CAL5</b>	66CAL5	250 ng /ml	2500 mIU/ml
1 ng/ml = 10 mIU/ml (1 <sup>st</sup> IRP 73/601)				

**Note:** The calibrators standards, human serum based, were calibrated using a reference preparation, which was assayed against the 1st International Reference Preparation (IRP# 73/601).

#### 4.3 CONJUGATE (code 66HRP)

CONJ|HRP

1 vial (13 mL) containing affinity purified antibody anti-CEA conjugated to HRP and biotinylated monoclonal IgG antibody (mouse) to CEA in buffer with dye and preservative.

Ready to use.

Store at 2-8°C.

#### 4.4 SUBSTRATE HS (LS.TMB-H2O2.HK)

SUBS|TMB

1 vial (13 mL) containing a stabilized mixture of TMB (3,3',5,5'-Tetramethylbenzidine) and H<sub>2</sub>O<sub>2</sub> (hydrogen peroxide). Cap label: blue.

Ready to use.

Store at 2-8°C.

#### 4.5 STOP SOLUTION (LS.STOP.K)

SOLN|STOP

1 vial (13 mL) of 0.3 M H<sub>2</sub>SO<sub>4</sub>. Cap label: red.

Ready to use.

Store at 2-8°C.

*The MSDS is available upon request of laboratory personnel.*

#### 4.6 WASHING BUFFER: (code TLAVM)(50X)

WASH|BUF|50X

1 bottle (20ml) of concentrated washing solution, containing a surfactant in phosphate buffered saline. A preservative has been added.

**Preparation of Wash Buffer working solution:** dilute contents of wash buffer concentrate 50x to 1000 mL with distilled or deionised water in a suitable storage container

**STORAGE** store at room temperature (20-27°C). Stable for up to 60 days.

#### 5.0 STORAGE AND STABILITY AFTER THE FIRST OPENING

- kit components should be stored at 2-8°C and should not be used after the expiry date.
- opened reagents are stable for 60 days when stored at 2-8°C.

#### 6.0 MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Distilled or deionized water and a clean one litre cylinder for the dilution of washing buffer.
- Pipette capable of delivering 25µl volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100ml and 0.300ml volumes with a precision of better than 1.5%.
- Microplate automatic washer or a squeeze bottle (optional).
- Microplate Reader with 450nm, 405 nm and 620nm wavelength OD capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate cover for incubation steps.
- Vacuum aspirator (optional) for wash steps.
- Timer.
- Quality control materials.

#### 7.0 WARNINGS OR PRECAUTIONS

##### 7.1 SAFETY PRECAUTION

- The kit is for in vitro diagnostic use only.
- Not for Internal or External Use in Humans or Animals
- All products that contain human serum have been found to be non reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA required tests. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease.
- Follow Good laboratory procedures for handling blood products.
- Avoid contact with reagents containing hydrogen peroxide, sulphuric and preservatives which may be toxic if ingested. Do not pipette by mouth.
- Follow all local, state and national regulations for disposal of all waste material.

##### 7.2 TECHNICAL PRECAUTION

- Avoid exposing TMB/H<sub>2</sub>O<sub>2</sub> reagent to direct sun-light, metal or oxidants.
- Avoid microbial contamination of reagents.
- The performance data represented here were obtained using specific reagents listed in the package insert.
- It is important that the time of reaction in each well is held constant for reproducible results.
- Pipetting of samples should not extend beyond 10 minutes to avoid assay drift.
- If more than 1 plate is used, it is recommended to repeat the dose response curve.
- Addition of the substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the addition of the substrate and the stopping solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.

- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.

#### 8.0 SPECIMEN COLLECTION AND STORAGE

- The specimens shall be blood, serum in type and the usual precautions in the collection of venipuncture samples should be observed.
- For accurate comparison to established normal values, a fasting morning serum sample should be obtained.
- The blood should be collected in a venipuncture tube without additives or anti-coagulants.
- Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.
- Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days.
- Avoid repetitive freezing and thawing.
- When assayed in duplicate, 0.050ml of the specimen is required.

#### 9.0 ASSAY PROCEDURE

##### 9.1 PREPARATION FOR ASSAY

- Bring all reagents to room temperature (20-27°C).
- Format the microplate wells for each calibrator, control and patient specimen to be assayed in duplicate. Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.

##### 9.2 PIPETTING AND INCUBATION STEPS

- A. Pipette** 25µl of the appropriate calibrator, control or specimen into the assigned wells
- B. Pipette** 100 µl of Conjugate into the wells;
- C. Swirl** the microplate gently for 20-30 seconds to mix and the **cover**.
- D. Incubate** at room temperature for **60 minutes**.
- E. Discard** the content of the microplate by aspiration or decantation. If decanting, blot the plate dry with absorbent paper;
- F. Wash** wells adding 300 µl of washing buffer working solution (see 4.6), then aspirate (or decant, tap and blot); Repeat 2 additional times for a total of 3 washes.
- An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage.*
- If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat 2 additional times.*
- G. Add** 100µl of Substrate HS to each well.
- H. Incubate** the wells at room temperature for **15 minutes**.
- I. Add** 100 µl of Stop Solution in the same sequence and timing adopted to dispense the Substrate HS
- J. Read** at 450 and 405 nm with the reference filter set at 620 nm within 30 minutes.

#### 10.0 CALCULATION OF RESULTS

##### 10.1 VALIDITY OF THE ASSAY

The OD of calibrator 5 should be  $\geq 1.3$ .

##### 10.2 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, medium and high ranges of the dose response curve for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

##### 10.3 OD CONVERSION

The optical densities (O.D.s) of some calibrators and samples may be higher than 2.0, in such a case, they could be out of the measurement range of the microplate reader. It is therefore necessary, for O.D.s higher than 2.0, to perform a reading at 405 nm (=wavelength of peak shoulder) in addition to 450 nm (peak wavelength) and 620 (reference filter for the subtraction of interferences due to the plastic).

For microplate readers unable to read the plate at 3 wavelengths at the same time,

it is advisable to proceed as follows:

- Read the microplate at 450 nm and at 620 nm.
- Read again the plate at 405 nm and 620 nm.
- Find out the wells whose ODs at 450 nm are higher than 2.0
- Select the corresponding ODs read at 405 nm and multiply these values at 405 nm by the conversion factor 3.0 (where OD 450/OD 405 = 3.0), that is:  $OD_{450\text{ nm}} = OD_{405\text{ nm}} \times 3.0$

Warning: The conversion factor 3.0 is suggested only. For better accuracy, the user is advised to calculate the conversion factor specific for his own reader.

**10.4 DATA REDUCTION - AUTOMATED METHOD**

Use the 4 parameters logistic – preferred – or the smoothed cubic spline function as calculation algorithm.

**NOTE:** If computer controlled data reduction is used to calculate the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.


**10.5 DATA REDUCTION: MANUAL METHOD**

A dose response curve is used to ascertain the concentration of CEA in unknown specimens.

- Record the OD obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the OD for each duplicate calibrator versus the corresponding CEA concentration in ng/ml on linear graph paper (do not average the duplicates of the calibrators before plotting).
- Draw the best-fit curve through the plotted points.
- To determine the concentration of CEA for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average OD (0.391) intersects the dose response curve at (22.5ng/ml) CEA concentration (See Figure 1).

EXAMPLE 1				
Sample ID	Well number	OD	Mean OD	Value ng/ml
CAL 0	A1	0.017	0.018	0
	B1	0.019		
CAL 1	C1	0.160	0.159	5
	D1	0.159		
CAL 2	E1	0.231	0.227	10
	F1	0.224		
CAL 3	G1	0.431	0.424	25
	H1	0.418		
CAL 4	A2	1.776	1.770	50
	B2	1.763		
CAL 5	C2	2.851	2.866	250
	D2	2.880		
patient	E2	0.389	0.391	22.5
	F2	0.384		

Figure 1

The data presented in Example 1 and Figure 1 are for illustration only and should not be used in lieu of a dose response curve prepared with each assay.

**10.6 MEASUREMENT RANGE**

Patient specimens with CEA concentrations above 250 ng/ml may be diluted (for example 1/10 or higher) with normal male serum (CEA < 5 ng/ml) and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor (10).

**11.0 EXPECTED VALUES**

Nearly 99% of non-smokers have CEA concentrations less than 5ng/ml. Similarly 99% of smokers have concentrations less than 10ng/ml (4).

TABLE I Expected Values for EIAGEN CEA kit	
Non-smokers	<5ng/ml
Smokers	<10ng/ml

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal"-persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

**12.0 LIMITATIONS OF THE PROCEDURE**

- Sample(s), which are contaminated microbiologically, should not be used in the assay. Highly lipemic or hemolysed specimen(s) should similarly not be used.
- CEA has a low clinical sensitivity and specificity as a tumor marker. Clinically an elevated CEA value alone is not of diagnostic value as a test for cancer and should only be used in conjunction with other clinical manifestations (observations) and diagnostic parameters. There are patients with colorectal cancer that do not exhibit elevated CEA values and elevated CEA values do not always change with progression or regression of disease. Smokers demonstrate a higher range of baseline values than non-smokers.

**13.0 PERFORMANCE CHARACTERISTICS**

**13.1 PRECISION**

The within and between assay precision of the EIAGEN CEA kit were determined by analyses on three different levels of control sera. The number, mean value, standard deviation (SD) and coefficient of variation for each of these control sera are presented in Table 3 and Table 4.

TABLE 2 Within Assay Precision (Values in ng/ml)				
Sample	N	MEAN	SD	C.V.
Level 1	20	4.8	0.35	7.3%
Level 2	20	21.7	1.35	6.2%
Level 3	20	60.5	3.58	5.9%

TABLE 3 Between Assay Precision* (Values in ng/ml)				
Sample	N	MEAN	SD	C.V.
Level 1	10	5.0	0.41	8.2%
Level 2	10	21.2	1.25	5.9%
Level 3	10	59.5	3.15	5.3%

\*As measured in ten experiments in duplicate

**13.2 ACCURACY-METHOD COMPARISON**

The EIAGEN CEA kit was compared with a reference Elisa method. Biological specimens from normal, and elevated concentrations were assayed. The total number of such specimens was 202. The least square regression equation and the correlation coefficient were computed for the EIAGEN CEA kit in comparison with the reference method. The data obtained is displayed in Table 4.

TABLE 4			
Method	Mean	Least Square Regression Analysis	Correlation Coefficient
EIAGEN CEA kit	5.67	y = -0.1164+1.0324(X)	0.998
Reference method	5.75		

**13.3 SENSITIVITY**

The EIAGEN CEA kit has a sensitivity of 0.025 ng. This is equivalent to a sample containing 1 ng/ml CEA concentration.

**13.4 SPECIFICITY**

Highly specific antibodies to CEA molecules have been used in the EIAGEN CEA kit. No interference was detected with the performance of EIAGEN CEA kit upon addition of massive amounts of the following substances to a human serum pool.

Acetylsalicylic Acid	100 µg/ml
Ascorbic Acid	100 µg/ml
Caffeine	100 µg/ml
AFP	10 µg/ml
PSA	1.0 µg/ml
CA-125	10,000 U/ml
hCG	1000 IU/ml
hLH	10 IU/ml
hTSH	100 mIU/ml
hPRL	100 µg/ml

**13.5 LINEARITY & HOOK EFFECT**

Three different lot preparations of the EIAGEN CEA kit reagents were used to assess the linearity and hook effect. Massive concentrations of CEA (> 60,000 ng/ml) were used for linear dilutions in pooled human patient sera. The test showed no hook effect up to concentrations of 60,000 ng/ml and a within dose recovery of 92.0 to 111.4%.

**14.0 AUTOMATION**

Application protocols for the proper automation on the Adaltis microtiterstrips analyzers (Labotech, Personal Lab and Nexgen) are available upon request at Adaltis directly.

## 15.0 SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors .

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different ( $\pm$ 50%) from OD reported on QC	<ul style="list-style-type: none"> <li>- incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes)</li> <li>- incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation)</li> <li>- error in washing or in photometer reading (suggestion: check operating or settings of respective instruments)</li> <li>- contamination of Substrate or Conjugate (suggestion: use only disposable and clean plastic containers)</li> </ul>
Low reproducible results	<ul style="list-style-type: none"> <li>- not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes)</li> <li>- error in washing or in reading (suggestion: check operating or settings of respective instruments)</li> <li>- contamination of Substrate (suggestion: use only disposable and clean plastic containers)</li> <li>- pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)</li> </ul>
no colourimetric reaction after addition of substrate	<ul style="list-style-type: none"> <li>- some reagent not pipetted</li> <li>- strong contamination of Conjugate or Substrate</li> <li>- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)</li> </ul>
too low reaction (too low ODs)	<ul style="list-style-type: none"> <li>- incubation time too short, incubation temperature too low</li> </ul>
too high reaction (too high ODs)	<ul style="list-style-type: none"> <li>- incubation time too long, incubation temperature too high</li> <li>- water quality for wash buffer insufficient (low grade of deionization)</li> <li>- insufficient washing (conjugates not properly removed)</li> </ul>
unexplainable outliers	<ul style="list-style-type: none"> <li>- contamination of pipettes, tips or containers</li> <li>- inconstant and insufficient washing (conjugates not properly removed)</li> </ul>
too high within-run CV%	<ul style="list-style-type: none"> <li>- reagents and/or strips not pre-warmed to Room Temp. prior to use</li> <li>- plate washer is not washing correctly (suggestion: clean washer head)</li> </ul>
too high between-run CV%	<ul style="list-style-type: none"> <li>- incubation conditions not constant (time, temperature)</li> <li>- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)</li> <li>- person-related variation</li> </ul>

## 15.0 REFERENCES

1. Gold P, Freedman SO, *J Exp Med*, **121**, 439 (1965).
2. Zamcheck N, *Adv Intern Med*, **19**, 413 (1974).
3. Rayncao G, Chu TM, *JAMA*, **220**, 381 (1972).
4. Wild D, *The Immunoassay Handbook*, Stockton Press (1994) p444.
5. Sorokin JJ, Sugarbaker PH, Zamcheck N, Pisick M, Kupchik HZ, Moore FD. "Serial carcinoembryonic antigen assays. Use in detection of cancer recurrence." *JAMA* 1974;228:49-53.
6. Mackay AM, Patel S, Carter S, Stecens U, Lawrence DJR, Cooper EH, et al. "Role of serial plasma assays in detection of recurrent and metastatic colorectal carcinomas." *Br. Med. J.* 1974; 4:382-385.
7. Sikorska H, Schuster J, Gold P. "Clinical applications of carcinoembryonic antigen." *Cancer Detection Preview* 1988; 12:321-355.
8. Minton JP, Martin EW Jr. "The use of serial CEA determinations to predict recurrence of colon cancer and when to do a secondlook surgery." *Cancer* 1978;42:1422-27.
9. Staab HJ, Anderer FA, Stumpf E, Fischer R. "Slope analysis of the postoperative CEA time course and its possible application as an aid in diagnosis of disease progression in gastrointestinal carcinoma." *Am. J.Surgery* 1978;136:322-327.
10. Thomas P, Toth CA, Saini KS, Jesup JM, Steele G. Jr., "The structure, metabolism and function of carcinoembryonic antigen gene family." *Biochem Biophys Acta* 1990;1032:177-189.
11. Yamashita K, Totami K, Kuroki M, Ueda I, Kobata A. "Structural studies of the carbohydrate moieties of carcinoembryonic antigens." *Cancer Research* 1987;47:3451-3459.
12. Hammerstrom S, Shively JE, Paxton RJ, Beatty BG, Larson A, Ghosh R, et al. "Antigenic sites in carcinoembryonic antigen." *Cancer Research* 1989;49:4852-58.
13. National Institute of Health. Carcinoembryonic Antigen: Its role as a marker in the management of cancer. A national Institute of Health Consensus Development Conference. *Ann.Intern.Med.*1981;94:407-409.

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
### 16.0 BIBLIOGRAFIA

- Gold P, Freedman SO, *J Exp Med*, **121**, 439 (1965).
- Zamcheck N, *Adv Intern Med*, **19**, 413 (1974).
- Raynao G, Chu TM, *JAMA*, **220**, 381 (1972).
- Wild D, *The Immunoassay Handbook*, Stockton Press (1994) p444.
- Sorokin JJ, Sugarbaker PH, Zamcheck N, Pisick M, Kupchik HZ, Moore FD. "Serial carcinoembryonic antigen assays. Use in detection of cancer recurrence." *JAMA* 1974;228:49-53.
- Mackay AM, Patel S, Carter S, Stecens U, Lawrence DJR, Cooper EH, et al. "Role of serial plasma assays in detection of recurrent and metastatic colorectal carcinomas." *Br. Med. J.* 1974; 4:382-385.
- Sikorska H, Schuster J, Gold P. "Clinical applications of carcinoembryonic antigen." *Cancer Detection Preview* 1988; 12:321-355.
- Minton JP, Martin EW Jr. "The use of serial CEA determinations to predict recurrence of colon cancer and when to do a secondlook surgery." *Cancer* 1978;42:1422-27.
- Staab HJ, Anderer FA, Stumpf E, Fischer R. "Slope analysis of the postoperative CEA time course and its possible application as an aid in diagnosis of disease progression in gastrointestinal carcinoma." *Am. J.Surgery* 1978;136:322-327.
- Thomas P, Toth CA, Saini KS, Jesup JM, Steele G. Jr., "The structure, metabolism and function of carcinoembryonic antigen gene family." *Biochem Biophys Acta* 1990;1032:177-189.
- Yamashita K, Totami K, Kuroki M, Ueda I, Kobata A. "Structural studies of the carbohydrate moieties of carcinoembryonic antigens." *Cancer Research* 1987;47:3451-3459.
- Hammerstrom S, Shively JE, Paxton RJ, Beatty BG, Larson A, Ghosh R, et al. "Antigenic sites in carcinoembryonic antigen." *Cancer Research* 1989;49:4852-58.
- National Institute of Health. Carcinoembryonic Antigen: Its role as a marker in the management of cancer. A national Institute of Health Consensus Development Conference. *Ann.Intern.Med.* 1981;94:407-409.

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







6. Mackay AM, Patel S, Carter S, Stecens U, Lawrence DJR, Cooper EH, et al. "Role of serial plasma assays indetection of recurrent and metastatic colorectal carcinomas." *Br. Med. Jr.* 1974; 4:382-385.
7. Sikorska H, Schuster J, Gold P. "Clinical applications of carcinoemryonic antigen." *Cancer Detection Preview* 1988; 12:321-355.
8. Minton JP, Martin EW Jr. "The use of serial CEA determinations to predict recurrence of colon cancer and when to do a secondlook surgery." *Cancer* 1978;42:1422-27.
9. Staab HJ, Anderer FA, Stumpf E, Fischer R. "Slope analysis of the postoperative CEA time course and its possible application as an aid in diagnosis of disease progression in gastrointestinal carcinoma." *Am. J.Surgery* 1978;136:322-327.
10. Thomas P, Toth CA, Saini KS, Jesup JM, Steele G. Jr. "The structure, metabolism and function of carcinoemryonic antigen gene family." *Biochem Biophys Acta* 1990;1032:177-189.
11. Yamashita K, Totami K, Kuroki M, Ueda I, Kobata A. "Structural studies of the carbohydrate moieties of carcinoemryonic antigens." *Cancer Research* 1987;47:3451-3459.
12. Hammerstrom S, Shively JE, Paxton RJ, Beatty BG, Larson A, Ghosh R, et al. "Antigenic sites in carcinoemryonic antigen." *Cancer Research* 1989;49:4852-58.
13. National Institute of Health. Carcinoemryonic Antigen: Its role as a marker in the management of cancer. A national Institute of Health Consensus Development Conference. *Ann.Intern.Med.* 1981;94:407-409.