EIAgen Ferritin KIT

REF LI4023

 96 tests

 **Adaltis Italia S.p.A.**

Via Cristoni, 12
40033 Casalecchio di Reno – (BO) Italy
Tel. +39-051-6136511 – Fax + 39-051-575280
www.adaltis.com

FOR IN VITRO DIAGNOSTIC USE ONLY

Store at 2...8 °C

en

SYMBOLS USED ON LABELS

MTPLATE	Microplate
CONJ HRP	Conjugate
CAL0	Calibrator 0
CAL1	Calibrator 1
CAL2	Calibrator 2
CAL3	Calibrator 3
CAL4	Calibrator 4
CAL5	Calibrator 5
SUBS TMB	Substrate HS (TMB)
SOLN STOP	Stop Solution
WASH BUF 20X	Wash buffer 20X
LOT	Lot number
REF	Catalogue code
	Expiry date (use by...)
 2°C  8°C	Temperature limitation (store at 2...8°C)
	Number of tests
	Keep away from sunlight
	Manufactured by
	Attention, See Instructions For Use
IVD	In vitro diagnostic medical device (In vitro diagnostic use)
	Biological risks (see section 7 WARNINGS AND PRECAUTIONS)

1.0 INTENDED USE

EIAgen Ferritin kit is a colourimetric immunoenzymatic method for the quantitative determination of the concentration of Ferritin in human serum and plasma.

2.0 SUMMARY AND EXPLANATION OF THE TEST

One of the most prevalent human disorders is the dietary deficiency of iron and the resulting anemia (1). The assays for iron, total iron-binding capacity, and other assessments of iron compounds in the body are clinically significant (2). Iron-storage compounds in the body include hemoglobin, hemosiderin, myoglobin, and the cytochromes. In most tissues, Ferritin is a major iron-storage protein (3). Human Ferritin has a molecular weight of approximately 450,000, and consists of a protein shell around an iron core. Each molecule of Ferritin may contain as many as 4,000 iron atoms (4). Under normal conditions, this may represent 25% of the total iron found in the body. In addition, Ferritin can be found in several isomers (5). High concentrations of Ferritin are found in the cytoplasm of the reticuloendothelial system, liver, spleen, and bone marrow (6). Methods previously used to measure iron in such tissues are invasive, cause patient trauma and lack adequate sensitivity (7). The measurement of Ferritin in serum is useful in determining changes in body iron storage, and is non-invasive with relatively little patient discomfort. Serum Ferritin levels can be measured routinely and are particularly useful in the early detection of iron-deficiency anemia in apparently healthy people (8, 9). Serum Ferritin measurements are also clinically significant in the monitoring of the iron status of pregnant women, blood donors, and renal dialysis patients. High Ferritin levels may indicate iron overload without apparent liver damage, as may be noted in the early stages of idiopathic hemochromatosis (10). Ferritin levels in serum have also been used to evaluate clinical conditions not related to iron storage, including inflammation, chronic liver disease, and malignancy (11).

3.0 PRINCIPLE OF THE ASSAY

EIAgen Ferritin kit is based on the simultaneous binding of human Ferritin to two monoclonal antibodies: one immobilized on microwell plates, the other conjugated with horseradish peroxidase (HRP).

After incubation the bound/free separation is performed by a simple solid-phase washing. The enzyme in the bound fraction reacts with the Substrate (TMB) developing a blue colour that turns into yellow after adding the Stop solution (H₂SO₄).

The Ferritin concentration in the sample is calculated on the basis of a series of calibrators.

The color intensity is proportional to the Ferritin concentration in the sample.

4.0 CONTENT OF THE KIT

4.1 MICROPLATE (code 15PIA)

MTPLATE

1 microplate, 12 x 8 breakable wells coated with monoclonal Anti-Ferritin antibodies.

Store at 2-8°C.

4.2 CALIBRATORS

6 vials containing human liver Ferritin with gentamicin 0.01% as preservative.

Open the package of coated microplate only when it is at room temperature and close immediately after use.

Store at 2-8°C.

Calibrator	Symbol	Code	Concentration*	Volume
Calibrator 0	CAL0	15CAL0	0 ng/mL	3 mL
Calibrator 1	CAL1	15CAL1	5 ng/mL	1.0 mL
Calibrator 2	CAL2	15CAL2	20 ng/mL	1.0 mL
Calibrator 3	CAL3	15CAL3	100 ng/mL	1.0 mL
Calibrator 4	CAL4	15CAL4	400 ng/mL	1.0 mL
Calibrator 5	CAL5	15CAL5	1000 ng/mL	1.0 mL

* approximate concentration.

The calibrators are calibrated against the WHO 1st IS Ferritin Human Liver 80/602

4.3 CONJUGATE (code 15HRP)

CONJ|HRP

1 vial containing 13 ml monoclonal anti-Ferritin antibody conjugate to horseradish peroxidase (HRP).

Ready to use.
Store at 2-8°C.

4.4 SUBSTRATE HS (LS.TMB-H202.HK)

SUBS|TMB

1 vial containing 13 ml of a stabilized mixture of TMB 0.25g/l (3,3',5,5'-Tetramethylbenzidine) and H₂O₂ (Hydrogen Peroxide).

Ready to use.
Store at 2-8°C.

4.5 STOP SOLUTION (LS.STOP.K)

SOLN|STOP

1 vial containing 13 ml of H₂SO₄ 0.3 M.

Ready to use.
Store at 2-8°C
The MSDS is available upon request of laboratory personnel

4.6 WASH BUFFER 20X (cod. TLAVD) **WASH|BUF|20X**

1 bottle containing 50 ml of concentrated washing buffer, containing phosphate buffer 50mM pH 7.4 and Tween20 1g/l.

Preparation of Wash Buffer working solution: dilute contents of wash buffer concentrate (50 ml) to 1000 ml with distilled water.

STORAGE store at room temperature (20-27°C). Stable until kit expiry date.

5.0 STORAGE AND STABILITY AFTER THE FIRST OPENING

- Kit components should be stored at 2-8°C and should be protected against sunlight.

6.0 MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Distilled water
- Automatic dispensers.
- Microplate reader capable of measuring the O.D.s at 450, 405 and 620 nm
- Precision pipettes: 0.02 and 0.1 ml
- Disposable pipette tips
- Absorbent paper
- Graph paper
- Control sera

7.0 WARNINGS AND PRECAUTIONS

7.1 SAFETY PRECAUTIONS

- Reagents contain Proclin 300^R and Gentamicin as preservatives.
- Avoid contact with reagents containing hydrogen peroxide, sulphuric acid and preservatives which may be toxic if ingested. Do not pipette by mouth.
- Calibrators contain human serum. They have been found to be non reactive for HbsAg and to anti-HIV antibodies by FDA; however, these products should be handled as potentially hazardous.
- Follow all local, state and national regulations for disposal of all waste material.

7.2 TECHNICAL PRECAUTIONS

- Follow the utmost precision when diluting and dispensing the reagents.
- Do not use reagents belonging to different lots.
- Reliable and reproducible results will be obtained when the assay procedure is carried out with a complete understanding of the package insert instructions and with adherence to good laboratory practice.
- The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.

8.0 SPECIMEN COLLECTION AND STORAGE

Ferritin determination can be carried out using either plasma or serum. If dispensing does not occur by five days from gathering, store specimen at -20°C.

9.0 ASSAY PROCEDURE

9.1 PREPARATION FOR ASSAY

- All reagents should be allowed to reach room temperature (18-25 °C) before use.
- All reagents should be mixed by gentle inversion or swirling prior to use. Do not induce foaming
- Since it is necessary to assay in duplicate, two wells should be prepared for each point of the Standard curve (0-5 Calibrators), two for each Specimen and one for the Blank.

9.2 PIPETTING AND INCUBATION STEPS

- Dispense** 20 µl of each calibrator and specimen in duplicate, leaving one well empty for the blank.
- Dispense** 100 µl of conjugate into the wells of calibrators and specimen; do not add conjugate into the Blank well.
- Incubate** 1 h at +22-28°C (room temperature).
- Remove** the reagent mixture, then **wash** dispensing 0.3 ml diluted wash solution into each well; wash twice again so as to completely remove the liquid.
- Dispense:** 100 µl Substrate into each well.
- Incubate** 10 minutes in the dark at room temperature (22+28°C).
- Dispense:** 100 µl Stop Solution into each well.
- Read** optical density, within 30 minutes, at 450 nm and 405 nm setting the reference filter at 620 nm. If it is possible to do it automatically, adjust the reader to zero using the Blank.

10.0 CALCULATION OF RESULTS

10.1 QUALITY CONTROL

Good laboratory practice requires that all controls (low, medium, and high) are run. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

10.2 OD CONVERSION

The optical densities (O.D.s) of some calibrators and samples may be higher than 2.0; in such a case, they could be out of the measurement range of the microplate reader. It is therefore necessary, for O.D.s higher than 2.0, to perform a reading at 405 nm (=wavelength of peak shoulder) in addition to 450

nm (peak wavelength) and 620 (reference filter for the subtraction of interferences due to the plastic).

For microplate readers unable to read the plate at 3 wavelengths at the same time, it is advisable to proceed as follows:

- Read the microplate at 450 nm and at 620 nm.
- Read again the plate at 405 nm and 620 nm.
- Find out the wells whose ODs at 450 nm are higher than 2.0
- Select the corresponding ODs read at 405 nm and multiply these values at 405 nm by the conversion factor 3.0 (where OD 450/OD 405 = 3.0), that is: $OD\ 450\ nm = OD\ 405\ nm \times 3.0$

Warning: The conversion factor 3.0 is suggested only. For better accuracy, the user is advised to calculate the conversion factor specific for his own reader.

10.3 DATA REDUCTION - AUTOMATED METHOD

Subtract the Blank OD from calibrator and sample ODs

Use the 4 parameters logistic – preferred – or the smoothed cubic spline function as calculation algorithm.

10.4 DATA REDUCTION - MANUAL METHOD

- Subtract the Blank OD from calibrator and sample ODs
- Calculate the average OD for each point in the standard curve and for each sample.
- Log-to-log draw on paper the values of Ferritin calibrators, then extrapolate the straight line.
- Read the corresponding Ferritin concentrations of controls and samples

10.5 MEASUREMENT RANGE

The present method allows the determination of Ferritin concentrations from 5 ng/ml to 1,000 ng/ml. Use the zero calibrator to dilute samples having concentrations > 1,000 ng/ml.

11.0 EXPECTED VALUES

Each laboratory should establish its own range on the basis of the patients' population.

The Ferritin plasma values are included in the following ranges:

Number of patients	Mean (ng/mL)	Range (ng/mL)
Fertile women:50	53	6 - 180
Women in post-menopause: 30	105	8 – 350
Men 50	175	20 - 400

12.0 LIMITATIONS OF THE PROCEDURE

- No "hook effect" was observed up to 50,000 ng/ml of Ferritin with the present method.
- No interferences from rheumatoid factors, bilirubin, hemoglobin, cholesterol, triglycerides have been observed with the present method.
- For diagnostic purposes, the results obtained by this test should be used as an adjunct to other data (e.g., symptoms, results of other tests, clinical impressions) available to the physician.

13.0 PERFORMANCE CHARACTERISTICS

13.1 PRECISION

a) Within-run precision.

Within-run precision has been evaluated by replicate determinations on three different control sera assayed in one assay. The within-run variability is reported below:

Serum sample	1	2	3
Number of replicates	16	16	16
Average Ferritin (ng/mL)	15.6	76.1	160.2
Standard Deviation	0.84	2.8	5.2
% CV	5.4	3.7	3.24

b) Between-run precision.

Between-run precision has been evaluated on three different control sera assayed in two different analyses. The between-run variability is reported below:

Serum sample	1	2	3
Number of replicates	16	16	16
Average Ferritin (ng/mL)	15.4	78.2	165.2
Standard Deviation	0.88	3.6	6.93
CV%	5.7	4.6	4.2

13.2 RECOVERY

Two patients were given known quantities of Ferritin in growing amount. The average recovery compared to the previous concentrations was 98.0 %.

	Expected Concentration (ng/ml)	Observed Concentration (ng/ml)	Recovery %
SAMPLE 1			
Q. added: 18			
50 ng/ml	68.0	65.0	95.6
100 ng/ml	118	122	103.4
200 ng/ml	218	210	96.3
SAMPLE 2			
Q. added: 128			
50 ng/ml	178	169	94.9
100 ng/ml	228	234	102.6
200 ng/ml	328	317	96.6
SAMPLE 3			
Q. added: 8.7			
12.5 ng/ml	21.2	22.4	106.2
25 ng/ml	33.7	31.8	94.4
50 ng/ml	58.7	56.4	96.1
100 ng/ml	108.7	105.9	97.4

13.3 LINEARITY

Two patient samples were serially diluted with zero calibrator for a linearity study. The mean recovery was 103.5 %.

	Expected Concentration (ng/mL)	Observed Concentration (ng/mL)	Expected%
Sample 1			
Dilution			
Non diluted		86	
1:2	43.0	45.1	105.4
1:4	21.5	22.1	102.8
1:8	10.75	12.4	107.9
Sample 2			
Dilution			
Non diluted		128	
1:2	64	65	101.6
1:4	32	30.4	93.7
1:8	16	17.1	106.1

13.4 METHOD COMPARISON

The ElAgen Ferritin kit (y) was compared to another RIA commercial method (x).

$$Y = 2.8 + 1.06X$$

$$r=0.99$$

13.5 SENSITIVITY

The minimal detectable concentration of Human Ferritin by this assay is estimated to be 1.0 ng/mL.

13.6 SPECIFICITY

The cross-reactivity of antibodies, as calculated at 50% according to Abraham, are shown in the following table:

Iso Liver Ferritin	100.0 %
Iso Spleen Ferritin	80.0 %
Iso Heart Ferritin	12.0 %

14.0 AUTOMATION

Application protocols for the proper automation on the Adaltis microtiterstrips analyzers (Labotech, Personal Lab and Nexgen) are available upon request at Adaltis directly.

15.0 SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors.

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different (\pm 50%) from OD reported on QC	<ul style="list-style-type: none"> - incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation) - error in washing or in photometer reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate or Conjugate (suggestion: use only disposable and clean plastic containers)
Low reproducible results	<ul style="list-style-type: none"> - not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - error in washing or in reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate (suggestion: use only disposable and clean plastic containers) - pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)
no colourimetric reaction	<ul style="list-style-type: none"> - some reagent not pipetted - strong contamination of Conjugate or Substrate

after addition of substrate	- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)
too low reaction (too low ODs)	- incubation time too short, incubation temperature too low
too high reaction (too high ODs)	<ul style="list-style-type: none"> - incubation time too long, incubation temperature too high - water quality for wash buffer insufficient (low grade of deionization) - insufficient washing (conjugates not properly removed)
unexplainable outliers	<ul style="list-style-type: none"> - contamination of pipettes, tips or containers - inconstant and insufficient washing (conjugates not properly removed)
too high within-run CV%	<ul style="list-style-type: none"> - reagents and/or strips not pre-warmed to Room Temp. prior to use - plate washer is not washing correctly (suggestion: clean washer head)
too high between-run CV%	<ul style="list-style-type: none"> - incubation conditions not constant (time, temperature) - controls and samples not dispensed at the same time (with the same intervals) (check pipetting order) - person-related variation

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17. Van Oost, S.A. Willekens, F.L.A., B.R. Van Neebors, B.R., Van Den Beld, B. Clin. Chem ,28/12 ,2429 - 2433 (1982).