EIAgen GI Marker CA 19-9 Kit

REF LI4040

 96 tests

 **Adaltis Italia S.p.A.**
Via Cristoni, 12
40033 Casalecchio di Reno – (BO) Italy
Tel. +39-051-6136511 – Fax +39-051-575280
www.adaltis.com

FOR IN VITRO DIAGNOSTIC USE ONLY

Store at 2...8 °C

en

SYMBOLS USED ON LABELS

MTPLATE	microplate
CONJ HRP 20x	HRP Conjugate Concentrate
BIOTIN Anti-CA19-9	Biotin Anti-CA 19-9
DIL CONJ	conjugate diluent
CAL0	Calibrator 0
CAL1	Calibrator 1
CAL2	Calibrator 2
CAL3	Calibrator 3
CAL4	Calibrator 4
CAL5	Calibrator 5
CONTROL1 CONTROL2	Controls
SOLN TMB	TMB substrate
SOLN STOP	Stop solution
WASH BUF 25X	washing buffer concentrate
LOT	Lot number
REF	Catalogue Code
	Expiry date (Use by...)
 2°C  8°C	Temperature limitation (store at 2...8°C)
	Number of tests
	Keep away from sunlight
	Manufactured by
	Attention, See Instructions For Use
IVD	In vitro diagnostic medical device (In vitro diagnostic use)
	Biological risks (see SECTION 7- Warnings and Precautions).

1.0 INTENDED USE

The EIAgen GI Marker CA19-9 kit is intended for the quantitative determination of the cancer associated antigen CA19-9 in human serum.

2.0 SUMMARY AND EXPLANATION OF THE TEST

The EIAgen GI Marker CA19-9 kit is based on a mouse monoclonal antibody, C192, highly specific for the sialyl Lewis^x epitope, also known as CA19-9 antigen (1). In adults, the epitope is typically expressed in trace amounts on mucosal cells of gastrointestinal epithelia.

In patients with malignant disease, the epitope may appear associated with high molecular weight mucin in blood. Assays for CA19-9 are frequently used to monitor gastrointestinal malignancies such as pancreatic, gall bladder, gastric, and colorectal cancers (2, 3).

3.0 PRINCIPLE OF THE TEST

The EIAgen GI Marker CA19-9 kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated together with a biotinylated Anti-CA19-9 monoclonal antibody (MAb) C192 (derived from mice) in streptavidin coated microtiter strips. CA19-9 present in calibrators or samples is adsorbed to the streptavidin coated microtiter strips by the biotinylated Anti-CA19-9 MAb during the incubation. The strips are then washed and incubated with HRP (Horseradish peroxidase) labeled Anti-CA19-9 MAb C192. After washing, buffered Substrate/Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetra-methylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue color will develop if antigen is present. The intensity of the color is proportional to the amount of CA19-9 present in the samples. The color intensity is determined in a microtiter plate photometer at 620 nm. Optionally the color turns to yellow after the addition of the Stop Solution and can be read at 405 nm.

Calibration curves are constructed for each assay by plotting optical density (OD) value versus the concentration for each calibrator. The CA19-9 concentrations of patient samples are then read from the calibration curve.

4.0 CONTENT OF THE KIT

4.1 MICROPLATE (code 35PIA) **MTPLATE**

1 microplate of 12 strips of 8 breakable The wells are coated with streptavidin. After opening, immediately return unused strips to the aluminium pouch, containing desiccant. Reseal carefully to keep dry.

4.2 CONJUGATE (code 35HRP) **CONJ|HRP|20X**

1 vial of 0.75 mL containing 20X concentrate Stock Solution of HRP conjugated to Anti-CA19-9 monoclonal antibody from mouse (approximately 40 µg/mL). Contains preservatives.

To be diluted with Conjugate Diluent prior to use.

Preparation of Conjugate (HRP-Anti-CA19-9) working solution: Prepare the required quantity of Conjugate working solution by mixing 50 µL of Conjugate 20x concentrate with 1 mL of Conjugate Diluent per strip (see table below).

N° of strips	CONJUGATE µl	Conjugate Diluent ml
1	50	1
2	100	2
3	150	3
4	200	4
5	250	5
6	300	6
7	350	7
8	400	8
9	450	9
10	500	10
11	550	11
12	600	12

To be sure of the stability use a clean plastic or glass bottle for preparation of the Conjugate working solution.

Alternative: Pour the content of the Conjugate 20x concentrate into the vial of Conjugate Diluent and mix gently. Make sure that the entire content of the Conjugate 20x concentrate is transferred to the vial of Conjugate Diluent.

STORAGE: the Conjugate working solution is stable for 3 weeks at 2–8°C in a sealed container. Do not prepare more Conjugate working solution than will be used within this period and make sure that it is stored properly.

4.3 CONJUGATE DILUENT (code 35HRPD) **DIL|CONJ**

1 vial containing 15 mL of phosphate citrate-buffer (pH 7.5) with bovine serum albumin, blocking agents, Tween 40, an inert blue dye, and 0.01 % methylisothiazolone (MIT) as preservative.

Ready for use.

4.4 CALIBRATORS

6 vials of 0.75 ml (except calibrator 0: 8ml). Each calibrator contains CA 19-9 antigen in a Tris-HCl buffered salt solution containing bovine serum albumin, an inert yellow dye and 0.05% sodium azide as preservative. Ready to use. The concentrations of calibrators are:

Calibrator	symbol	code	concentration
Calibrator 0	CAL0	35CAL0	0 U/mL
Calibrator 1	CAL1	35CAL1	15 U/mL
Calibrator 2	CAL2	35CAL2	30 U/mL
Calibrator 3	CAL3	35CAL3	60 U/mL
Calibrator 4	CAL4	35CAL4	120 U/mL
Calibrator 5	CAL5	35CAL5	240 U/mL

Calibrator 0 should also be used for dilution of samples as Sample Diluent.

Reference materials.

Since no common reference material is available for CA19-9 antigen, ADALIS CA19-9 EIA Calibrator values are assigned against a set of in-house reference standards.

4.5 BIOTIN Anti CA 19-9 (code 35BT) BIOTIN Anti-CA19-9

1 vial containing 15 ml of solution of biotin labelled Anti-CA19-9 monoclonal antibody from mouse (approximately 2 µg/mL). Contains phosphate citrate buffer (pH 7.5), bovine serum albumin, blocking agents, Tween 40, an inert red dye, and 0.01 % methylisothiazolone (MIT) as preservative.

Ready to use.

4.6 CONTROLS

Control	symbol	code	
Control 1	CONTROL 1	35CTRL1	1 vial of 0.75 mL
Control 2	CONTROL 2	35CTRL2	1 vial of 0.75 mL

1 vial of each control containing 0.75 mL of CA 19-9 antigen in a tris-HCl buffered salt solution containing bovine serum albumin and 0.05% sodium azide as preservative.

Ready to use.

4.7 TMB-SUBSTRATE (code TMBB)

SOLN TMB

1 vial containing 12 mL of ready-to-use Hydrogen peroxide and 3,3',5,5' tetra-methylbenzidine (TMB).

Ready to use.

4.8 STOP SOLUTION (code STOPB)

SOLN STOP

1 vial containing 15 mL of 0.12 M hydrochloric acid.

Ready to use.

4.9 WASHING BUFFER 25X (code TLAVB)

WASH BUF 25X

1 bottle containing 50 mL of wash buffer concentrate. Wash buffer concentrate contains a Tris-HCl buffered salt solution with Tween 20. contains German II as preservative.

Preparation of Wash Buffer working solution: Pour the 50 mL Washing Buffer Concentrate into a clean container and dilute 25-fold by adding 1200 mL of distilled or deionised water to give a Wash Buffer working solution.

STORAGE: the stability of diluted wash buffer is 2 weeks at 2-25°C in a sealed container.

5.0 STORAGE AND STABILITY AFTER THE FIRST OPENING

- Store kit components at 2...8° C and do not use after the reported expiry date.
- After use, the microplate should be resealed, the bottle caps replaced and tightened and the kit stored at 2...8°C.

The opened kit should be used within expiry date indicated on labels.

6.0 MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Microtiter plate shaker:** Shaking should be medium to vigorous. Longitudinal shaking approximately 200 strokes/min, oscillations 700—900/min.
- Microtiter plate wash device:** Automatic plate wash capable of performing 1, 3 and 6 washing cycles, or semi manual microtiter plate washing device connected to vacuum pump or water-jet vacuum and a liquid trap for retaining aspirated liquid. The Nunc Immuno-8 manual strip washer is recommended if an automatic microtiter plate wash is not used.
- Microtiter plate photometer:** with a wavelength of 620 nm and/or 405 nm, and an OD range of 0 to 3.0.
- Precision pipettes:** with disposable plastic tips for dispensing microlitre volumes. An 8-channel pipette or respenser pipette with disposable plastic

tips for delivery of 100 µL is useful but not essential. Pipettes for dispensing millilitre volumes.

- Distilled or deionised water:** For preparation of Wash Solution.
- Others:** reagent reservoirs for multichannel pipettes, paper towels, and timer

7.0 WARNINGS OR PRECAUTIONS

7.1 SAFETY PRECAUTIONS

- All reagents in this kit are for in vitro diagnostic use only
- Only experienced laboratory personnel should use this test and handling should be agreement with GLP.
- Operators should wear gloves and protective clothing when handling any patient sera or serum based products.
- Avoid contact with reagents containing hydrogen peroxide, hydrochloric acid and preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is total of all detergents.
- Materials used in the preparation of human source reagents have been tested for the presence of antibodies to Human Immunodeficiency Virus (HIV 1 and 2), as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.
- Follow all local, state and national regulations for disposal of all waste material. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.

7.2 TECHNICAL PRECAUTIONS

A. Correct use of reagents.

- The performance data represented here were obtained using specific reagents listed in this package insert. Do not use reagents from other manufacturers in the kits.
- Do not use reagents from other EIAgen kits with this kit. Do not mix reagents from different kit lots.
- Do not dilute or adulterate the kit reagents, unless directed by the kit protocol.
- The **TMB Substrate** is very sensitive for contamination. The TMB HRP Substrate should be colourless or slightly bluish. A blue colour indicates that the reagent has been contaminated and should be discarded. For optimal stability of the TMB Substrate, pour the required amount from the vial to a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or respenser pipette tip).

B. Correct pipetting procedure.

- Be sure to use clean disposable plastic pipette tips and a proper pipetting technique when handling samples and reagents. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid. A proper pipetting technique is of particular importance when handling the TMB HRP-Substrate solution.

C. Correct washing procedure.

- Pre-washing of the strips before the assay is necessary
 - Microplate washing is important. Improperly washed wells will give erroneous results
- WASHING PROCEDURE:** A careful washing procedure of the strips is essential. Ensure that each well is filled up completely to the top edge and that the aspiration of the wells between and after the washing cycles is complete and the wells are dry. If there is liquid left in the wells, invert the plate and tap it carefully against absorbing paper.
- Automatic strip washer:** Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles prior to and after each incubation step. The aspiration/wash device should not be left standing with the Wash Solution for long periods as the needles may get clogged, giving poor liquid delivery and suction.

D.. Adherence to assay procedure and specifications

- The test protocol must be followed strictly. Observe the indicated incubation times and temperature and the washing procedure. These are critical steps. Do not allow the wells to dry between incubations. Please refer also to section 9.0_ASSAY PROCEDURE.
- Include the positive and negative controls in every test run to monitor for reagent stability and correct assay performance. The obtained OD values for calibrators and concentration for controls have to be always compared to the ones reported in QC sheet. Please refer also to section 10.1_VALIDITY OF THE ASSAY AND 10.2_QUALITY CONTROL.
- Do not use the kit to determine values outside the range indicated in the Instruction for Use. Please refer also to section 10.5_MEASUREMENTS RANGE.

8.0 SPECIMEN COLLECTION AND STORAGE

The kit is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Samples can be stored at 2–8°C for 24 hours. For longer periods stores samples at -70°C or below. Samples should not be stored in a self-defrosting freezer and not be thawed and refrozen before analysis. Allow frozen samples to thaw slowly, preferably at 2–8°C over night and then bring the samples to room temperature before analysis.

9.0 ASSAY PROCEDURE

9.1 REAGENT PREPARATION

- Bring all reagents to room temperature (20...25°C). The assay should only be performed at temperatures between 20–25°C to obtain accurate results.
- Prepare washing buffer working solution according to par. 4.9;
- Prepare conjugate working solution according to par 4.2;
- Select sufficient microwells for the test. Transfer the required number of microtiter plate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant and reseal carefully).
- Before starting to pipette calibrators and patient specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
- Perform each determination in duplicate for both calibrators and patient samples. A calibration curve should be run with each assay.

9.2 PIPETTING AND INCUBATION STEPS

- A. Wash** each strip once with the Wash Buffer working solution. Do not wash more strips than can be handled within 30 min.
- B. Pipette** in duplicate 25 µL of the Calibrators (CAL 0, 1, 2, 3, 4, 5), Controls (C 1, C 2) and patient specimens into the strip wells;
- C. Add** 100 µL of Biotin Anti-CA19-9 to each well using a 100 µL precision pipette (or an 8-channel 100 µL precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
- D. Incubate** the plate for 2 hours (± 5 min) at room temperature (20-25°C) with constant shaking of the plate using a microtiter plate shaker.
- E.** After the first incubation **aspirate** and **wash** each strip 3 times using the wash procedure described in section 7.2C
- F. Add** 100 µL of Conjugate working solution to each well. Use the same pipetting procedure as in item E above.
- G. Incubate** the frame for 1 hour (± 5 min) at room temperature (20–25°C) with constant shaking.
- H.** After the second incubation **aspirate** and **wash** each strip 6 times, using the wash procedure described in section 7.2C.
- I. Add** 100 µL of TMB HRP-Substrate to each well using the same pipetting technique as in item 4. The TMB Substrate should be added to the wells as quickly as possible and the time between addition to the first and last well should not exceed 5 min.
- J. Incubate** for 30 min (± 5 min) at room temperature with constant shaking. Avoid exposure to direct sunlight.
- K.** Immediately **read** the OD at 620 nm in a microtiter plate photometer.

Option If the laboratory does not have access to a microtiter plate reader capable of reading at 620 nm, the OD can be determined in this way: Add 100 µL of Stop Solution, mix and read the OD at 405 nm in a microtiter plate photometer within 15 min after addition of Stop Solution.

10.0 CALCULATION OF RESULTS

10.1 VALIDITY OF THE ASSAY

Control 1 and 2 may be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated.

10.2 QUALITY CONTROL

Each laboratory may in addition prepare its own serum pools at different levels, which can be used as internal controls in order to assure the precision of the assay.

If the results should fall out with the expected range, repeat the assay using freshly prepared controls. If the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory

procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall out with the acceptable ranges.

10.3 DATA REDUCTION: AUTOMATED METHOD

If a microtiter plate photometer with built-in data calculation program is used, refer to the manual for the photometer and create a program using the concentration stated on the label of each of the CA19-9 calibrators.

For automatic calculation of CA19-9 results it is recommended to use either of the following methods:

- **Cubic spline** curve fit method. Calibrator 0 should be included in the curve with the value 0 U/mL.
- **Spline smoothed curve** fit method. Calibrator 0 should be used as plate blank.
- **Interpolation with point-to-point evaluation.** Calibrator 0 should be included in the curve with the value 0 U/mL.
- **Quadratic curve** fit method. Calibrator 0 should be included in the curve with the value 0 U/mL.

NOTE: 4-parametric or linear regression should not be used.

10.4 DATA REDUCTION: MANUAL METHOD

A calibration curve is constructed by plotting the OD values obtained for each Calibrator against the corresponding CA19-9 concentration (in U/mL), see figure below. The unknown CA19-9 concentrations can then be read from the calibration curve using the mean OD value of each patient specimen.

10.5 MEASUREMENT RANGE

The EIAgen Gi MARKER CA19-9 kit measures concentrations between 1 and 240 U/mL. If CA19-9 concentrations above the measuring range are to be expected, it is recommended to dilute samples with the Calibrator 0 prior to analysis.

If samples in an initial analysis give CA19-9 levels higher than 240 U/mL the samples should be diluted 1/10 and 1/100 with CA19-9 Calibrator 0 to obtain the accurate CA19-9 concentration of the samples.

1/10 dilution = 50 µL of specimen + 450 µL of CA19-9 Calibrator 0

1/100 dilution = 50 µL of 1:10 dilution + 450 µL of CA19-9 Calibrator 0

The CA19-9 concentration of the undiluted sample is then calculated as:

Dilution 1/10= 10 x measured value

Dilution 1/100= 100 x measured value

11.0 EXPECTED VALUES

CA19-9 was measured in 100 healthy blood donors, 36 women and 64 men. The mean value obtained was 6.5 U/mL with a standard deviation of 6.4. The median value was 4.5 U/mL, range 0 – 29 U/mL. The lower and upper extremes of the normal range were examined using IFCC recommended non-parametric statistical treatment. The reference interval contains the central 95% fraction of the reference distribution. Reference limits may accordingly be estimated as the 2.5% (lower) and 97.5% (upper) fractiles. These limits cut off a fraction of 2.5% of the values in each tail of the reference distribution. Non-parametric estimates: N=100

Fraction	Reference limit (U/mL)
2.5th (lower)	0
97.5th (upper)	25

100% of the healthy subjects had assay values below 37 U/mL. It is recommended that each laboratory establish its own normal range to account for such local environmental factors as diet, climate, living conditions, patient selection, etc.

12.0 LIMITATION OF THE PROCEDURE.

The level of CA19-9 cannot be used as absolute evidence for the presence or absence of malignant disease and the CA19-9 test should not be used in cancer screening. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the management of patients, and the CA19-9 test should not replace any established clinical examination.

Benign conditions such as acute or chronic pancreatitis or cholelithiasis may cause elevated levels of CA19-9.

Patients with a Le a-/b- phenotype do not express the CA19-9 reactive epitope. Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffers.

13.0 PERFORMANCE CHARACTERISTICS

13.1 PRECISION

Total precision was determined according to NCCLS guideline EP5-A (4) using four levels of frozen pooled human serum containing added ascitespool. Each sample was randomly pipetted (n=2/analysis) and analysed twice each day

over 20 days. The analyses were undertaken during a period of 53 months, by > three different technicians and using 20 different kit batches.

Sample	Replicates	Mean U/mL	Within-run SD (U/mL)	Within run CV %	Between day SD (U/mL)	Between day CV %
CA19-9 1	80	15.4	0.6	3.8	1.0	6.8
CA19-9 2	80	56.3	1.9	3.3	3.6	6.3
CA19-9 3	80	99.8	4.5	4.5	6.2	6.2
CA19-9 4	80	182	7.9	4.4	12	7.0

13.2 DETECTION LIMIT

The detection limit of the ADALTIS CA19-9 EIA assay is < 1 U/mL defined as the concentration corresponding to the mean of the OD values for the Calibrator 0 plus 2 standard deviations according to the formula:

$$(2 \times \text{SD CAL 0}) / (\text{OD CAL 1} - \text{OD CAL 0}) \times 15 \text{ U/ml}$$

13.3 RECOVERY

Spiked serum samples were prepared by adding different levels of human CA19-9 antigen to normal serum samples. The recovery of the antigen was in the range 90–110%.

13.4 HOOK EFFECT

No hook effect has been noticed with samples up to 50 000 U/mL.

NOTE: In very high samples the colour of the substrate will change from blue to greenish (and eventually yellow in extremely high samples). This will lead to a falsely low OD at 620 nm, and in extreme cases the OD may fall within the calibration curve range and noticed as a hook.

13.5 LINEARITY

Eight patient samples were serially diluted with Calibrator 0 and analyzed. The obtained values were between 96-105 % of the expected values in the range 10–200 U/mL.

13.6 SPECIFICITY

The monoclonal antibody used (C192) is highly specific for the sialyl Lewis^a epitope (1). The NCCLS guideline EP7-P (5) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

	Concentration with no significant (± 10%) interference
Lipemia (Intralipid®)	10 mg/mL
Bilirubin, unconjugated	0.6 mg/mL
Hemoglobin	5 mg/mL

14.0 AUTOMATION

Application protocols for the proper automation on the Adaltis microtiterstrips analyzers are available upon request at Adaltis directly.

15.0 SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors .

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different (± 50%) from OD reported on QC	<ul style="list-style-type: none"> - incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation) - error in washing or in photometer reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate (suggestion: use only disposable and clean plastic containers)
Low reproducible results	<ul style="list-style-type: none"> - not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - error in washing or in reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate (suggestion: use only disposable and clean plastic containers) - pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)
no colourimetric reaction after addition of substrate	<ul style="list-style-type: none"> - some reagent not pipetted strong contamination of conjugate or Substrate - errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)
too low reaction (too low ODs)	<ul style="list-style-type: none"> - incorrect conjugate dilution - incubation time too short, incubation temperature too low
too high reaction (too high ODs)	<ul style="list-style-type: none"> - incorrect conjugate dilution - incubation time too long, incubation temperature too high - water quality for wash buffer insufficient (low grade of deionization) - insufficient washing (conjugates not properly removed)
unexplainable outliers	<ul style="list-style-type: none"> - contamination of pipettes, tips or containers - inconstant and insufficient washing (conjugates not properly removed)
too high within-run CV%	<ul style="list-style-type: none"> - reagents and/or strips not pre-warmed to Room Temp. prior to use - plate washer is not washing correctly (suggestion: clean washer head)
too high between-run CV%	<ul style="list-style-type: none"> - incubation conditions not constant (time, temperature) - controls and samples not dispensed at the same time (with the same intervals) (check pipetting order) - person-related variation

16.0 REFERENCES

1. Rye P. D. et al., (1998) Summary report on the ISOBM TD-6 workshop: Analysis of 20 monoclonal antibodies against Sialyl Lewis^a and related antigens. Montreaux, Switzerland; September 19-24 1997. *Tumor Biol.* 19: 390-420.
2. Eskelinen M., and Haglund U. (1999) Developments in serologic detection of human pancreatic adenocarcinoma. *Scand J Gastroenterol.* 34: 833-844.
3. Hammarström S., and Stigbrand T. (2002) Gastric Cancer. In "Tumor markers, Physiology, Pathobiology, Technology and Clinical Applications" Eds. Diamandis E. P. et al., AACCC Press, Washington pp 375-379.
4. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999)
5. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986)

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16.0 RIFERIMENTI

- Rye P. D. et al., (1998) Summary report on the ISOBM TD-6 workshop: Analysis of 20 monoclonal antibodies against Sialyl Lewis x and related antigens. Montreux, Switzerland; September 19-24 1997. *Tumor Biol.* 19: 390-420.
- Eskelinen M., and Haglund U. (1999) Developments in serologic detection of human pancreatic adenocarcinoma. *Scand J Gastroenterol.* 34: 833-844.
- Hammarström S., and Stigbrand T. (2002) Gastric Cancer. In "Tumor markers, Physiology, Pathobiology, Technology and Clinical Applications" Eds. Diamandis E. P. et al., AACC Press, Washington pp 375-379.
- National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999)
- National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986)

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16.0 BIBLIOGRAFIA

1. Rye P. D. et al., (1998) Summary report on the ISOBM TD-6 workshop: Analysis of 20 monoclonal antibodies against Sialyl Lewis_x and related antigens. Montreaux, Switzerland; September 19-24 1997. *Tumor Biol.* 19: 390-420.
2. Eskelinen M., and Haglund U. (1999) Developments in serologic detection of human pancreatic adenocarcinoma. *Scand J Gastroenterol.* 34: 833-844.
3. Hammarström S., and Stigbrand T. (2002) Gastric Cancer. In "Tumor markers, Physiology, Pathobiology, Technology and Clinical Applications" Eds. Diamandis E. P. et al., AACC Press, Washington pp 375-379.
4. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999)
5. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986)