ImmuGlo™ AUTOANTIBODY TEST SYSTEM

For *in vitro* Diagnostic Use

| | | | |
|------------|-------|---------|---|
| REF | Code: | 1100 | Mouse Liver Substrate 48 Determinations |
| REF | Code: | 1107 | COMVI Mouse Kidney/Stomach Substrate 48 Determinations |
| REF | Code: | 1136 | COMVI-III Mouse Kidney/Stomach/Liver Substrate 48 Determinations |
| REF | Code: | 1136-96 | COMVI-III Mouse Kidney/Stomach/Liver Substrate 96 Determinations |

PRODUCT INSERT

INTENDED USE

Indirect immunofluorescence (IF) antibody tests for the detection and quantitation of anti-nuclear antibodies (ANA) *Code: 1100, 1107, 1136, 1136-96*, anti-mitochondrial antibodies (AMA), anti-smooth muscle antibodies (ASMA) and anti-gastric parietal cell antibodies (AGPA) *Catalog No. 1107, 1136, 1136-96*, in human serum.

SUMMARY AND EXPLANATION

Antinuclear antibodies (ANA), detected by indirect immunofluorescence, aid in the diagnosis of connective tissue disorders including systemic lupus erythematosus (SLE), mixed connective tissue disease, Sjögren's syndrome and scleroderma¹⁻⁵. ANA occur in about 95% of SLE patients as well as patients with other connective tissue diseases. ANA may also occur in other disorders such as chronic active hepatitis and primary biliary cirrhosis⁶⁻⁸.

Anti-mitochondrial antibodies (AMA) occur in over 90% of primary biliary cirrhosis cases, 3-11% of chronic active hepatitis patients and are absent in patients with extra-hepatic biliary obstruction and in other liver diseases. The universal presence of anti-mitochondrial antibodies in primary biliary cirrhosis and their virtual absence in extra-hepatic jaundice makes their detection of considerable value in the differential diagnosis⁶⁻¹².

Anti-smooth muscle antibodies (ASMA) in high titer (≥ 160) occur in the majority of cases of chronic active hepatitis and in intermediate titers (40-80) in acute viral hepatitis. Occasionally they may occur in cases of primary biliary cirrhosis where they are also found in intermediate titers. The significance of titers of 20-40 is doubtful since these titers may occur in normal individuals^{13,14}.

Anti-gastric parietal cell antibodies (AGPA) are commonly associated with pernicious anemia and chronic atrophic gastritis where they occur in about 90% and 50% of cases, respectively. However, they are not disease specific as they may occur in low frequency in other disorders. Although healthy individuals may have gastric parietal cell antibodies, this finding may reflect asymptomatic atrophic gastritis. Negative findings for gastric parietal cell antibodies provide strong evidence for excluding pernicious anemia¹⁵⁻¹⁷.

PRINCIPLES OF PROCEDURE

In the indirect immunofluorescence method used in this kit, patients' sera are incubated on mouse kidney/stomach sections to allow binding of antibodies to the substrate. Any antibodies not bound are removed by rinsing the slide. Bound antibodies of the IgG class are detected by incubation of the substrate with fluorescein-labeled, anti-human IgG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of ANA, ASMA, AMA and AGPA is demonstrated by an apple green fluorescence of specific histologic structures in the tissue. The titer (the reciprocal of the highest dilution giving a positive reaction) is then determined by testing serial dilutions ¹⁸.

PRODUCT INFORMATION

Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

Materials provided

| | | |
|------------|--|-------------------|
| REF | 1100 Mouse Liver Substrate | 48 Determinations |
| REF | 1107 COMVI Mouse Kidney/Stomach Substrate | 48 Determinations |
| REF | 1136 COMVI-III Mouse Kidney/Stomach/Liver Substrate | 48 Determinations |
| REF | 1136-96 COMVI-III Mouse Kidney/Stomach/Liver Substrate | 96 Determinations |

| | | |
|------------|--|--|
| 6 | SORB SLD 8 | 8 well Substrate Slides (1100, 1107, 1136) |
| 12 | SORB SLD 8 | 8 well Substrate Slides (1136-96) |
| 1 x 0.5 ml | CONTROL + ANA * | ANA Positive Control. Contains human serum. |
| 1 x 0.5 ml | CONTROL + AMA * | AMA Positive Control. Contains human serum. |
| 1 x 0.5 ml | CONTROL - * | Negative Control. Contains human serum. |
| 1 x 5 ml | IgG-CONJ FITC * | Anti-human IgG FITC Conjugate. Protect from light. (1136-96 contains 2 x 5ml) |
| 1 x 5 ml | IgG-CONJ FITC EB ** | Anti-human IgG FITC Conjugate containing Evan's Blue. Protect from light. (1136-96EB contains 2 x 5ml) |
| 1 x 60 ml | BUF * | Buffered Diluent. |
| 2 vials | BUF WASH | Phosphate Buffered Saline (PBS). Dissolve each vial to 1 liter. |
| 1 x 5.0 ml | MOUNTING MEDIUM * | Mounting Medium. Do not freeze. |
| 1 x 1.0 ml | EVANS | Evan's Blue Counterstain. |
| 1 x 12 | COVER SLD | Coverslips. |

* Contains < 0.1% NaN₃

† Replaces conjugate without counterstain in Code numbers containing "EB"

Material required but not provided

Fluorescence microscope
Micropipette or Pasteur pipette
Serological pipettes
Staining dish (e.g. Coplin jar)
Small test tubes (e.g. 13 x 75 mm) and test tube rack
Distilled or deionized water
1 liter container
Wash bottle
Paper towels
Incubation chamber

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials¹⁴.

WARNING - Sodium azide (NaN_3) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

SPECIMEN COLLECTION AND PREPARATION

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

PROCEDURE

Test Method

A. Screening

1. Dilute each patient serum 1:10 (Code 1107, 1136, 1136-96) with the Buffered Diluent provided (20 μl serum + 180 μl Diluent) or 1:20 (Code 1100) (10 μl serum + 190 μl Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 μl of the Negative Control to well #1. Similarly apply 1 drop of ANA Positive Control to well #2. If applicable, 1 drop of AMA Positive Control to well #3 (Code: 1107, 1136, 1136-96). Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 μl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.

7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 μ l) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evans blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:10 (Code: 1107, 1136, 1136-96) or 1:40 (Code: 1100). The reciprocal of the highest dilution producing a positive reaction is the titer.

Preparation of Serial Dilutions

Number four tubes 1 through 6. Add 0.9 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

| Tubes | 1 | 2 | 3 | 4 | 5 | 6 |
|------------------|--------|--------|--------|--------|--------|------------|
| Serum | 0.1 ml | | | | | |
| | + | | | | | |
| Buffered Diluent | 0.9 ml | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml |
| | | ↻ | ↻ | ↻ | ↻ | ↻ |
| Transfer | | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml |
| Final dilution | 1:10 | 1:20 | 1:40 | 1:80 | 1:160 | 1:320 etc. |

QUALITY CONTROL

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence of the nuclei, smooth muscle, tubules of the

kidney or gastric parietal cells. The AMA Positive Control should have 2+ or greater staining intensity of the tubules of the kidney. The ANA Positive Control should have 2+ or greater staining intensity of the nuclei of the kidney and liver with a predominantly homogeneous pattern.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

INTERPRETATION OF RESULTS

The results of the tests for ANA, AMA, ASMA, and AGPA antibodies should be reported as negative (< 10) on kidney and stomach sections (*Code: 1107, 1136, 1136-96*), negative (<20) on liver sections (*Code: 1100*), or alternatively, positive with titer.

Read only fields which contain specific staining of the nuclei for ANA, the kidney tubules for AMA, the blood vessel walls for ASMA and gastric parietal cells only for AGPA. All other reactions should be reported as negative for ANA, AMA, ASMA and/or AGPA.

ANA can be detected on all substrates but should be quantified on the kidney or HEp-2 cells. The nuclear staining patterns observable with the kidney substrate or HEp-2 cells provided include homogeneous, peripheral (rim), speckled and nucleolar. The centromere staining pattern (including mitotic figures) is seen most easily on HEp-2 cells. These nuclear staining patterns are described below. They may be one or a combination of several staining patterns. The latter are due to reactions to several different nuclear antigens.

- Homogeneous:** The entire nucleus fluoresces evenly with a diffuse staining pattern.
- Peripheral (rim):** The nuclear membrane stains most intensely with decreasing staining intensity of the nucleoplasm towards the center of the nucleus.
- Speckled:** Discrete coarse to fine round speckles fluoresce throughout the nucleus.
- Nucleolar:** The nucleoli stain as multiple solid bodies within the nucleus.

The specificity of some of the antibodies giving the above staining patterns may be further identified by tests for antibodies to nDNA and to various extractable nuclear antigens. These may be of diagnostic significance as listed in FIG 1 at the end of this document.

AMA may be observed on both the distal and proximal tubules of the kidney with the distal tubules staining more brightly. Even though the cytoplasm of the gastric parietal cells also stains, AMA should be quantitated on the kidney.

Staining of the stomach muscularis and kidney glomeruli may also be observed with ASMA, but only ASMA seen on the blood vessel walls of the kidney should be reported.

LIMITATION OF THE PROCEDURE

In some cases, sera positive for ANA may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titers determined.

In some cases the presence of two or more antibodies in a serum which are reactive with

the same substrate may cause an interference in their detection by immunofluorescence. This interference may cause either failure to detect ANA or suppression of its titer if the interfering antibody has a higher titer than ANA. All ANA reactions should be reported.

EXPECTED VALUES

As seen in Tables 1, 2, 3, 4 and 5 at the end of this document, tests for nuclear antibodies are used to screen for SLE and certain other immunologic disturbances. AMA occur in over 90% of cases of primary biliary cirrhosis and 3-11% of cases of chronic hepatitis. ASMA occur in the majority of cases of chronic active hepatitis and AGPA are commonly associated with pernicious anemia and chronic atrophic gastritis.

PERFORMANCE CHARACTERISTICS

The ImmuGlo™ Autoantibody Test System (Mouse Kidney/Stomach Sections) was compared with another commercially available fluorescent antibody test using mouse kidney/stomach as a substrate. The comparison included: 20 samples of ANA positive sera, 19 samples of AMA positive sera, 19 samples of ASMA positive sera, 20 samples of AGPA positive sera and 38 serum samples from normal subjects. Sera were tested starting at a 1:10 dilution with the procedure recommended by the manufacturer. These yielded comparable results as summarized in Tables 6 and 7 at the end of this document.

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REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI

1. Tan EM. Autoantibodies to nuclear antigens (ANA): Their immunobiology and medicine. *Adv Immunol* 33: 167-240, 1982.
2. Kumar V, Beutner EH and Chorzelski TP. Autoimmunity of the skin. In "Concepts in Immunopathology", Vol 1, Cruse JM and Lewis RE Jr, Eds, Karger, Basel, 318-353, 1985.
3. Reimer G, Cornell RC and Tan EM. The biochemical nature of nuclear antigens reactive with antinuclear antibodies. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 519-531, 1987.
4. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Standardization of antinuclear antibody and other immunofluorescent tests used in immunopathologic studies of the skin. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 41-64, 1987.
5. Tan EM, Chan EKL, Sullivan KF and Rubin RL. Antinuclear antibodies (ANAs): Diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity. *Clin Immunol Immunopathol* 47: 121-141, 1988.
6. Manns M, Gerken G, Kyriatsoulis A and Meyer zum Büschenfelde KH. Significant autoimmune markers of autoimmune liver disorders: Current status. *J Clin Lab Anal* 1: 362-370, 1987.
7. Mackay IR. Autoimmunity and the liver. *Clin Aspects Immunity* 2: 8- 17, 1988.
8. McMillan SA, Alderdice JM, McKee CM et al. Diversity of autoantibodies in patients with anti-mitochondrial antibody and their diagnostic value. *J Clin Path* 4: 232-236, 1987.
9. Gershwin ME, Coppel RL and Mackay IR. Primary biliary cirrhosis and mitochondrial autoantigens - insights from molecular biology. *Hepatology* 8: 147-151, 1988.
10. Berg PA and Klein R. Mitochondrial antigens and autoantibodies from anti-M1 to anti-M9. *Klin Wochenschr* 64: 897-909, 1986.
11. Popper H and Paronetto F. Clinical, histologic and immunopathologic features of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 339-354, 1980.
12. Berg PA and Bacon H. Serology of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 355-373, 1980.
13. Anderson P, Small JV and Sobieszek A. Studies on the specificity of smooth muscle antibodies. *Clin Exp Immunol* 22: 22-29, 1975.
14. Kurki P, Miettinen A, Linder E, Pikkarainen P, Vuorio M and Salaspuro MP. Different types of smooth muscle antibodies in chronic active hepatitis and primary biliary cirrhosis: Their diagnostic and prognostic significance. *Gut* 21: 878-884, 1980.
15. Fisher JB and Taylor KB. The significance of gastric antibodies. *Brit J Haematol* 20: 1-7, 1971.
16. Chisholm M. Immunology of gastritis. *Clin Gastroenterol* 5: 419-428, 1976.
17. Bigazzi PE, Burek CL and Rose NR. Antibodies to tissue-specific endocrine, gastrointestinal and neurological antigens. In "Manual of Clinical Laboratory Immunology". Rose NR, Friedman H and Fahey JL, Eds, American Society for Microbiology, Washington DC, 762-770, 1986.
18. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
19. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].
20. Nisengard RJ. Antinuclear antibodies: Significance of titers. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Bean S, Eds, John Wiley and Sons, New York, 2nd Ed, 387-398, 1979.
21. Meyer zum Büschenfelde KH, Manns M and Trautman F. Autoimmunity in chronic liver diseases - relationship to SLE? In "Recent Advances in Systemic Lupus Erythematosus". Lambert PH, Perrin L, and Izui S, Academic Press, New York, 259-269, 1984.
22. Walker JG, Doniach D, Roitt IM and Sherlock S. Serologic tests in the diagnosis of primary biliary cirrhosis. *Lancet* 1: 827, 1965.
23. Paronetto F and Popper H. Hetero-iso- and autoimmune phenomena in the liver. In "Textbook of Immunopathology", Miescher PA and Müller-Eberhard HJ, Eds, Grune and Stratton, New York, 2nd Ed, 789-817, 1976.

Figure 1. Diagnostic Significance of Antinuclear Antibodies

| IF Staining Pattern | Nature of Antigen | Associated Disease |
|----------------------------|------------------------------|----------------------------|
| Homogeneous | Deoxyribonucleoprotein | SLE with renal involvement |
| Peripheral | DNA | SLE |
| Speckled | RNP | SLE or MCTD* |
| | Sm | SLE |
| | SS-A/SS-B | SLE or Sjögren's |
| | Scl-70 | Scleroderma |
| | 4S-6S RNA probably U3 RNA | Scleroderma |

*Mixed Connective Tissue Disease

Table 1: Incidence of Antinuclear Antibodies (ANA) Detected by Indirect Immunofluorescence on HEp-2 Cells

| Clinical Condition | No. of Sera | % Positive |
|------------------------------|--------------------|-------------------|
| SLE | 12 | 100 |
| Subacute Cutaneous LE (SCLE) | 7 | 86 |
| Scleroderma | 6 | 100 |
| Rheumatoid Arthritis | 10 | 50 |
| Normal Controls | 15 | 0 |

Table 2: Incidence of Antinuclear Antibodies (ANA) Detected by Indirect Immunofluorescence on Mouse Kidney Substrate

| Clinical Condition | No. of Sera | % Positive |
|---------------------------|--------------------|-------------------|
| SLE | 21 | 95 |
| Scleroderma | 17 | 82 |
| Rheumatoid Arthritis | 20 | 5 |
| Normal Controls | 96 | 0 |

Table 3: Incidence of Anti-Mitochondrial Antibodies (AMA) Detected by IFA on Mouse Kidney Substrate

| Clinical Condition | % Positive |
|--|-------------------|
| Primary Biliary Cirrhosis | 100 |
| Autoimmune Chronic Active Hepatitis | 8 |
| HBsAg and Chronic Active Hepatitis | 0 |
| Extrahepatic Jaundice and Other Liver Diseases | 0 |
| Systemic Lupus Erythematosus (SLE) | 3 |
| Rheumatoid Arthritis | 0 |
| Normal Controls | 0 |

Adapted from Meyer zum Büschenfelde KH, et al.²¹; Walker JG, et al.²² and Paronetto F and Popper H²³.

Table 4: Incidence of Anti-Smooth Muscle Antibodies (ASMA) as Detected by IFA on Mouse Kidney Substrate

| Clinical Condition | % Positive |
|-----------------------------------|-------------------|
| Chronic Active Hepatitis (Type A) | 50-87 |
| Primary Biliary Cirrhosis | 25 |
| Acute Viral Hepatitis | 87 |
| Infectious Mononucleosis | 87 |
| Burkitt's Lymphoma | 73 |
| Nasopharyngeal Carcinoma | 23 |
| Hodgkin's Disease | 23 |
| Myeloproliferative Disorder | 5 |
| Warts | 4 |
| Normal Controls | 3-18 |

Adapted from Anderson P, et al.¹³

Table 5. Incidence of Anti-Gastric Parietal Cell Antibodies (AGPA) as Detected by Indirect Immunofluorescence on Mouse Stomach Substrate

| Clinical Condition | % Positive |
|---------------------------------------|------------|
| Pernicious Anemia (PA) | 85-95 |
| Chronic Atrophic Gastritis without PA | 30-60 |
| Gastric Ulcer | 25-30 |
| Autoimmune Endocrinopathies | 25-33 |
| Sjögren's Syndrome | 30 |
| First Degree Relatives of PA Patients | 30 |
| Normal Controls | |
| < 20 years old | 2 |
| 20-60 years old | 6-8 |
| > 60 years old | 16 |

Table 6. Findings in Positive Sera

| | n | Negatives | | Positives | | |
|---------------------------|----|-----------|-------|----------------|---------|----------|
| | | < 10 | 10-20 | Titer 40-80 | 160-320 | 640-2560 |
| ANA Positive Sera | | | | | | |
| IMMCO | 20 | 0 | 0 | 7 | 6 | 7 |
| Other | 20 | 0 | 0 | 7 | 8 | 5 |
| AMA Positive Sera | | | | | | |
| IMMCO | 19 | 4 | 1 | 3 | 1 | 10 |
| Other | 19 | 4 | 1 | 3 | 4 | 7 |
| ASMA Positive Sera | | | | | | |
| IMMCO | 19 | 3 | 4 | 8 | 4 | 0 |
| Other | 19 | 2 | 5 | 6 | 5 | 1 |
| AGPA Positive Sera | | | | | | |
| IMMCO | 20 | 0 | 1 | 6 | 4 | 9 |
| Other | 20 | 0 | 2 | 6 | 7 | 5 |

Table 7. Findings in Normal Controls

| | n | Negatives | | Positives | |
|---------------------------|----|-----------|-------|----------------|------|
| | | < 10 | 10-20 | Titer 40-80 | >160 |
| ANA Positive Sera | | | | | |
| IMMCO | 38 | 30 | 3 | 5 | 0 |
| Other | 38 | 36 | 0 | 2 | 0 |
| AMA Positive Sera | | | | | |
| IMMCO | 38 | 38 | 0 | 0 | 0 |
| Other | 38 | 38 | 0 | 0 | 0 |
| ASMA Positive Sera | | | | | |
| IMMCO | 38 | 35 | 2 | 1 | 0 |
| Other | 38 | 29 | 8 | 1 | 5 |
| AGPA Positive Sera | | | | | |
| IMMCO | 38 | 37 | 0 | 1 | 0 |
| Other | 38 | 37 | 0 | 1 | 0 |

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
 Autorizzato/Representante Autorizado/Représentant Autorisé
 EMERGO Group, Inc.
 Molenstraat 15, 2513 BH, The Hague,
 The Netherlands
 Tel (+31) 345 8570, Fax (+31) 346 7299
 www.emergogroup.com

For technical assistance please contact:

IMMCO Diagnostics, Inc.

60 Pineview Drive

Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

IMMCO
DIAGNOSTICS

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immcodiagnostics.com

or your local product distributor