# ImmuGlo™ Anti-Keratin Antibody IFA

**REF** 1122 48 Determinations

## INTENDED USE

An indirect immunofluorescence antibody test for the detection and quantitation of anti-keratin antibodies (AKA) in human serum.

## SUMMARY AND EXPLANATION

Rheumatoid arthritis (RA) is the most prevalent systemic rheumatic disease afflicting 1-2% of the population. It is characterized by mononuclear cell infiltration and proliferation of synovial cells. The disease results in membrane inflammation followed by irreversible degradation of joint cartilage and bone structure. The pathogenesis of RA is unknown, however, it is characterized by the presence of various circulating autoantibodies such as to rheumatoid factor (RF), anti-keratin antibodies (AKA), and anti-perinuclear factor (APF). RF is present in 70-90% of patients with RA and is included in the ARA\* classification criteria.

Anti-keratin antibodies, initially described by Young et al<sup>1</sup>, are found to be highly specific for RA. AKA can be detected by indirect immunofluorescence (IFA) on rat esophagus substrate, even prior to the onset of joint symptoms<sup>2-6</sup>. They occur in approximately 40% of patients with RA, of which approximately 14% are RF negative. The association between RF and AKA is close. Circulating immune complexes are found in significantly higher concentrations in RA patients positive for AKA. This may explain the association of AKA with severe forms of RA.

## PRINCIPLES OF PROCEDURE

In the indirect immunofluorescence method used in this kit, patients' sera are incubated on rat esophagus sections to allow binding of antibodies to the tissue substrate. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgG class are detected by incubation of the substrate with fluorescein-labeled, anti-human IgG. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of AKA demonstrated by an apple green fluorescence of the stratum corneum of the mucosal epithelium. The titers (the reciprocal of the highest dilution giving a positive reaction) are then determined by testing serial dilutions.

## PRODUCT INFORMATION

### Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

### Materials provided

IMMCO Anti-Keratin Antibody IFA **REF** 1122

Kits contain sufficient reagents to perform 48 determinations each.

6 x	<b>SORB</b> <b>SLD</b> <b>8</b>
1 x 0.5 ml	<b>CONTROL</b> <b>+</b> <b>AKA</b> *
1 x 0.5 ml	<b>CONTROL</b> <b>-</b> *

### 8 well Rat Esophagus Slides

**Positive Control.** Contains human serum with BSA.

**Negative Control.** Contains human serum with BSA.

1 x 5 ml **IgG-CONJ** **FITC** \*

1 x 60 ml **BUF** \*

2 vials **BUF** **WASH**

1 x 5.0 ml **MOUNTING MEDIUM** \*

1 x 1.0 ml **EVANS**

1 x 12 **COVER** **SLD**

**Anti-human IgG FITC Conjugate.**  
Contains BSA. **Protect from light.**

**Sample Diluent.** Contains BSA.

**Phosphate Buffered Saline (PBS).**  
Dissolve each vial to 1 liter.

**Mounting Medium.** Do not freeze.

**Evan's Blue Counterstain.**

**Coverslips.** Reconstitute to one liter each.

\* Contains < 0.1% NaN<sub>3</sub>

### Material required but not provided

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes
- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation chamber

### WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>14</sup>.

WARNING - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

### SPECIMEN COLLECTION AND PREPARATION

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

### PROCEDURE

#### Test Method

#### A. Screening

1. Dilute each patient serum 1:10 with the Sample Diluent provided (0.1 ml serum + 0.9 ml Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.

4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 0.05 ml) of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 0.05 ml) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 0.05 ml) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evans blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

## B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:10. The reciprocal of the highest dilution producing a positive reaction is the titer.

### Preparation of Serial Dilutions

Number 6 tubes 1 through 6. Add 0.9 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4
Serum	0.1 ml			
	+			
Buffered Diluent	0.9 ml	0.2 ml	0.2 ml	0.2 ml
		↗	↗	↗
Transfer		0.2 ml	0.2 ml	0.2 ml
Final dilution	1:10	1:20	1:40	1:80 etc.

## **QUALITY CONTROL**

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence, whereas the Positive Control should have 2+ or greater staining intensity.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include: improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

## **INTERPRETATION OF RESULTS**

The results of the tests for AKA should be considered negative (< 10), positive (greater or equal to 320), or alternatively, positive with titer.

Read for staining of the stratum corneum layer of the epithelium for AKA, as shown in figure 1 at the end of this document.

## **LIMITATIONS OF THE PROCEDURE**

In some cases, sera positive for AKA may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titers determined.

Two or more types of autoantibodies in the serum may react with heart substrate causing interference with their detection and correct identification by immunofluorescence. This results either in the failure to detect AKA antibodies or suppression of their titer, if the interfering antibody has a higher titer.

## **EXPECTED VALUES**

AKA are found in sera from patients with pre-RA and RA. They are rarely detected in non-RA patients.

Specificity:	99%
Sensitivity:	44%
Positive Predictive Value:	43%

The incidence of AKA in relation to RA appears in Table 1 at the end of this document.

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
## REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI

1. Young BJJ, Mallya RK, Leslie RDG, Clark CJM, Hamblin TJ. Antikeratin antibodies in rheumatoid arthritis. *Br Med J*; 1979, ii:97-99.
2. Aho K, von Essen R, Kurki P, Paluso T, Heliövaara M. Antikeratin antibody and antiperinuclear factor as markers for subclinical rheumatoid disease process. *J Rheumatol*; 1993, 20:1278-1281.
3. Kurki P, Aho K, Paluso T, Heliövaara M. Immunopathology of rheumatoid arthritis. Antikeratin antibodies precede the clinical disease. *Arthr Rheum*; 1992, 35:914-917.
4. Paimela L, Gripenberg M, Kurki P, Leirisalo-Repo M. Antikeratin antibodies: diagnostic and prognostic markers for early rheumatoid arthritis. *Ann Rheum Dis*; 1992, 51:743-746.
5. von Essen R, Kurki P, Isomäki H, Okubo S, Kautiainen H, Aho K. Prospect for an additional laboratory criterion for rheumatoid arthritis. *Scand J Rheumatol*; 1993, 22:267-272.
6. Vincent C, Serre G, Lapeyre F, Fournié B, Ayrolles C, Fournié A, Soleilhavoup J. High diagnostic value in rheumatoid arthritis of antibodies to the stratum corneum of rat oesophagus epithelium, so-called 'antikeratin antibodies'. *Ann Rheum Diseases*; 1989, 48:712-722.
7. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].

**Figure 1: AKA Positive Reaction**

**Table 1: AKA and RF in Patients With Peripheral Oligo/Polyarthritis in Relation to ARA\* Criteria for RA**

AKA	≥ 4 Criteria		3 Criteria		< 3 Criteria	
	RF+	RF-	RF+	RF-	RF+	RF-
<b>Positive</b>	27%	4%	18%	9%	4%	0%
<b>Negative</b>	44%	23%	18%	54%	12%	84%

Scand J Rheumatol; 1993, 22:267-272.

\*American Rheumatology Association

*For technical assistance please contact:*

**IMMCO**<sup>®</sup>  
DIAGNOSTICS

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive**

**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)