# OTOblot™ Anti-68 kD (hsp-70)\* Western Blot

For Research Use only

## PRODUCT INSERT

**Code: 1190**      **20 Determinations**

### INTENDED USE

A Western Blot assay for the detection and identification of antibodies to the 68 kD (hsp-70) antigen associated with sensorineural hearing loss (SNHL).

### SUMMARY AND EXPLANATION

Hearing loss may be caused by a number of conditions. Certain types of hearing loss can be reversed if diagnosed early and appropriate treatment is instituted. Sensorineural hearing loss (SNHL), commonly referred to as nerve deafness, may be caused by genetic or acquired factors such as infections or can be immunologically mediated. Correct diagnosis of SNHL can be made with the aid of a combination of comprehensive patient history and laboratory studies. In the majority of cases, no cause of SNHL is apparent: such cases are referred to as idiopathic SNHL. A subgroup of idiopathic SNHL cases are treatable with immunosuppressive therapy with gratifying results<sup>1,2</sup>. The laboratory work up to identify these cases should include serum antibody tests to 68 kD (hsp-70) inner ear antigen<sup>3-8</sup>. Twenty-two percent of patients with bilateral rapidly progressive SNHL and 30% of the patients with Meniere's disease had antibodies that reacted with 68kD antigen present in the bovine inner ear extracts<sup>7</sup>. Anti-68 kD (hsp-70) antibodies also occur in approximately 60% of patients with bilateral and 35% of patients with unilateral Meniere's Syndrome and 37% of patients with contralateral delayed endolymphatic hydrops. In situations where corticosteroids are contraindicated, methotrexate or cytoxan treatment can be prescribed<sup>9</sup>.

### PRINCIPLE OF PROCEDURE

This Western Blot immunoassay uses purified recombinant inducible hsp-70 antigen (rhsp-70) from the bovine kidney. During the manufacturing process, recombinant purified hsp-70 kD antigen is mixed with a molecular weight alignment marker of 116 kD. This mixture is then subjected to polyacrylamide gel electrophoresis and transferred to a PVDF membrane. These membranes are cut into 3 mm x 5 cm strips for immunostaining.

To perform the test, strips are incubated with diluted patient serum. Antibodies specifically bind to the rhsp-70 protein on the strip. After proper washing and two incubation steps using Conjugate A and Conjugate B, strips are washed and incubated with a precipitable substrate. Anti-hsp-70 antibody positive reactions appear as a blue-violet band at 70kD.

### REAGENTS

#### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use, if liquid reagents are turbid or a precipitate is present. Prior to starting the assay, reagents must be equilibrated to room temperature (~22°). Antigen strips can only be used once. Do not interchange components of different lots. Do not use reagents beyond expiration date indicated on labels.

#### Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by US FDA required tests. However, all human blood derivatives and patient specimens should be considered potentially infectious and good laboratory practices in storing, dispensing and disposing of these materials must be followed<sup>10</sup>.

WARNING - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. NaN<sub>3</sub> is toxic if ingested. Report incidents immediately to laboratory director or poison control center.  
Follow good laboratory practices to minimize microbial and cross contamination of reagents.

**Materials provided**

**Code: 1190**

The OTOblot anti-68 kD (hsp-70) antibody Western Blot kit\*, contains sufficient reagents to perform 20 determinations.

- 1 x 20 Western Blot Strips, **containing rhsp-70.**
- 1 x 110 µl \*Positive Control (**red vial cap**)
- 1 x 110 µl \*Negative Control (**yellow vial cap**).
- 1 Control Card
- 1 x 240 µl \*Conjugate A (**blue vial cap**)
- 1 x 240 µl \*Conjugate B (**white vial cap**)
- 1 x 60 ml \*Serum Diluent
- 1 x 25 ml \*Enzyme Substrate (**amber bottle**)
- 1 vial \*Powdered Wash Buffer; **reconstitute to one liter with deionized or distilled water.**
- 3 Assay trays
- 2 Report Forms

**\*Contains <0.1% NaN<sub>3</sub>**

**Materials Required But Not Provided**

- Clean 1000 ml graduated cylinder
- Non-serrated forceps (Filter forceps)
- Rocker or rotating platform shaker
- Absorbent paper or paper towels
- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 10 to 1000 µl
- Disposable pipet tips
- Timer

**SPECIMEN COLLECTION AND HANDLING**

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

**PROCEDURE**

**Procedural Notes**

- Read instructions carefully before starting with the assay.
- Let serum specimens and test reagents equilibrate to room temperature for ~30 minutes prior to starting the test procedure. Return all unused specimens and reagents to the refrigerator promptly after use.
- Proper washing technique is critical to the satisfactory performance of the assay.
- Manipulate test strips with clean forceps only. Do not touch with bare hands.
- Strips are individually numbered at the top of each strip. Assign specimen identification numbers to the respective strips on the Report Form.
- Complete all other relevant information on the Report Form prior to starting the assay.

## Test Method

1. Using blunt forceps, place required number of **Strips** labeled side up into individual wells of the assay tray.
2. Pipet **1.0 ml** of Serum Diluent into each well.
3. Pipet **10 µl** of Positive and Negative Control and patient sample into appropriate wells to obtain a **1:101 dilution**. Incubate **60 minutes** ( $\pm 5$  min.) at room temperature on a rocker or rotating shaker.
4. Aspirate sample solution into waste container. Thoroughly wash strips with Wash Buffer by squirting approximately 2ml of solution directly onto strips. Wash strips with gentle agitation for **5 minutes** and aspirate solution into waste container. **Repeat 3x**. *Caution: Complete washing of the strips between incubations is crucial to obtain valid results. Improper washing will result in high background staining.*
5. Pipet **1.0 ml** of Serum Diluent followed by **10 µl** of Conjugate A into each well. Incubate **30 minutes** ( $\pm 5$  min) at room temperature on rocker or rotating shaker.
6. Repeat **Step 4**.
7. Pipet **1.0 ml** of Serum Diluent followed by **10 µl** of Conjugate B into each well. Incubate **30 minutes** ( $\pm 5$  min.) at room temperature on rocker or rotating shaker.
8. Repeat **Step 4**.
9. Pipet **1.0 ml** Substrate into each well and incubate with gentle shaking **20 minutes** ( $\pm 5$  min) at room temperature and reduced light.
10. Repeat **Step 4**, washing twice instead of four times.
11. Using blunt forceps, remove strips from assay tray and place them gently onto absorbent paper. Handle strips only at the ends and let them dry **15-20 minutes**.

## Quality Control

Though the control cards are lot specific, negative and positive controls must be included in each test run to ensure proper performance of the assay.

*Positive Control Reaction:* The positive control strip will exhibit a dense blue-violet band at 70kD (anti-rhsp-70 band) and the molecular weight alignment marker at 116 kD (**Figure 1**).

*Negative Control Reaction:* Typically, the negative control will exhibit only one band at the 116 kD molecular weight alignment marker. Any weak band at the 70 kD level on the negative control strip should be **considered the baseline reaction** for that particular test run (**Figure 2**).

## RESULTS

### Reading and Interpretation Guidelines

The OTOblot™ strips contain two proteins, one at 116 kD and the other is rHSP70 at 70 kD (**Figure 1**). The 116 kD protein serves as a molecular weight alignment marker.

1. Hold test strip between the positive and negative control reactions on the provided, laminated Control Card and align test strip using the 116 kD molecular weight alignment marker as the reference point (**Figure 2**). A 116 kD marker line on the control cards facilitates accurate alignment of the test strips with the control strips.
2. Compare reaction of test strip with those of the controls on either side. Use of a magnifying glass can facilitate proper alignment of test strips along the molecular weight marker and assist in observation of weak reactions.
3. Check for a crisp band on the test strip that aligns with the 70 kD band on the positive control strip. Such a reaction should be considered positive. Positive reactions can occur in varying intensities from weak to strong. Weak reactions should be compared with baseline reaction intensities at the corresponding position on the negative control strip. Reactions that are distinctly darker or denser than the baseline intensity at 70kD should be considered positive. Occasionally, minor lower molecular weight bands may be visible with certain sera. However, such reactions are not relevant and should be considered as negative for antibodies to 68kD (rhsp-70) antigen.

**Note:** Automated Western blot machines can occasionally show extra bands that may not be related to the 70kD band. Ensure the washing cycle is thorough. Repeat process manually for clarity.

## EXPECTED VALUES

Antibodies to 68 kD (hsp-70) occur in patients with active idiopathic SNHL and antibody titers are shown to correlate with disease activity. See tables 1 and 2 at the end of this document.

Studies evaluating 34 patients with rapidly progressive hearing loss demonstrated that anti-rhsp-70 OTOBlot had 42% sensitivity, 90% specificity, and 91% positive predictive value for predicting steroid responsiveness<sup>12</sup>.

## LIMITATIONS OF THE PROCEDURE

The OTOblot anti-68 kD (hsp-70) antibody Western Blot assay should be used as an aid to diagnosis. Positive results may be found in other autoimmune conditions or certain infectious diseases. Hence results should be evaluated and interpreted by a medical authority in light of the patient's clinical history and other laboratory findings. Some sera may react to the MW marker occasionally, the significance of which is not known.

## TROUBLESHOOTING GUIDE

- **Strong band at 68-70kD on Negative Control strip.** Likely cause: contaminated Negative Control vial, or cross contamination from well containing a positive serum.
- **Positive Control appears like Negative Control strip.** Likely cause: Negative Control vial was confused as Positive Control vial.
- **Strips are completely blank.** Likely cause: addition of Conjugate (A or B) or Substrate was omitted.
- **High background and poor contrast between bands and background.** Likely cause: wash step(s) may have been omitted or incorrectly performed, or incubations were overextended.

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

## REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI

1. Harris EN, Phil M, Asherson RA et al. Antiphospholipid antibodies - antibodies with a difference. *Ann Rev Med*; 1988, 39:261-271.
2. Harris EN, Gharavi AE, Wasley GD and Hughes GRV. Use of an enzyme-linked immunosorbent assay and inhibition studies to distinguish between antibodies to cardiolipin from patients with syphilis or autoimmune disorders. *J Infect Dis*; 1988, 157:23-31.
3. Lockshin MD. Anti-cardiolipin antibody. *Arth Rheum*; 1987, 30:471-472.
4. Derksen RHW, Beisma D, Bouma BN et al. Discordant effects of prednisone on anticardiolipin antibodies and the lupus anti-coagulant. *Arth Rheum*; 1986, 29:1295-1296.
5. Alarcón-Segovia D, Deleze M, Oria CV et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. *Medicine*; 1989, 68:353-365.
6. Triplett DA. Antiphospholipid antibodies and recurrent pregnancy loss. *Am J Reprod Immunol*; 1989, 20:52-67.
7. Montecucco C, Longhi M, Caporali R et al. Hematological abnormalities associated with anticardiolipin antibodies. *Hematologica*; 1989, 74:195-204.
8. Asherson RA, Khamashta MA, ORDI-Ros J, et al. The "primary" antiphospholipid syndrome: major clinical and serological features. *Medicine*; 1989, 68:366-374.
9. Walport MJ. Pregnancy and antibodies to phospholipids. *Ann Rheum Dis*; 1989, 48:795-797.
10. Machworth-Young CG, Loizou S and Walport MJ. Antiphospholipid antibodies and disease. *Medicine*; 1989, 72:767-777.
11. Harris EN, Gharavi AE, Hughes GRV. Antiphospholipid antibodies. *Clin Rheumatol*; 1985, 1:591-609.
12. Lubbe WF, Palmer SJ, Butler WS and Laggins GC. Fetal survival after prednisone suppression of maternal lupus anticoagulant. *Lancet*; 1983, 1:1361-1363.
13. Hogan MJ, Brunet DG, Ford PM et al. Lupus anticoagulant, antiphospholipid antibodies and migraine. *Can J Neurological Sci*; 1988, 15:420-425.
14. Harris EN, Pierangeli S. Antiphospholipid antibodies: method of detection. *Am J Reproduct Immunol*; 1992, 28:208-210.
15. Silveira LH, Lopez LR, Uzzle RT et al. IgA, IgG and IgM anticardiolipin antibodies in black patients with systemic lupus erythematosus. *Arthr Rheum 56th Annual Meeting*; 1992, S125.
16. Harris, NE. The second international anti-cardiolipin standardization workshop/The Kingston antiphospholipid antibody study (KAPS) group. *Am J Clin Path*; 1990, 94:476-484.
17. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health; 1993, (HHS Pub. No [CDC] 93-8395).
18. Alarcón-Segovia D, Delezé M, Oria CV et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. A perspective analysis of 500 consecutive patients. *Medicine*; 1987, 68:353-365.

**Figure 1**

**Figure 2**

**Table 1. Anti-68 kD (hsp-70) Antibodies in Patients with Idiopathic Bilateral SNHL (11)**

Disease Group	No. Tested	No. Positive	% Positive
IPBSNHL	72	42	58
Otosclerosis	11	0	0
Cogan's Syndrome	8	0	0
Normals	53	1	2

**Table 2. Correlation of Anti-68 kD (hsp-70) Antibody Reactivity to Disease Activity (11)**

Disease Activity	Anti-68 kD Antibody
Active	89%
Inactive	0%

**NOTES:**

*For technical assistance please contact:*

**OTOIMMUNE Diagnostics™**  
*a division of IMMCO Diagnostics, Inc.*  
**60 Pineview Drive**  
**Buffalo, NY 14228-2120**  
**Telephone: (716) 691-0091**  
**Fax: (716) 691-0466**  
**Toll Free USA/Canada: 1-800-537-TEST**  
**E-mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.  
Molenstraat 15, 2513 BH, The Hague,  
The Netherlands

REV.DEC2003

Tel (+31) 345 8570, Fax (+31) 346 7299  
[www.emergogroup.com](http://www.emergogroup.com)

Document No. PI4190