EN

ImmuliTM Enhanced dsDNA Antibody ELISA

Double Stranded DNA Antibody ELISA

IVD For *in vitro* diagnostic use

PRODUCT INSERT

REF 5120 dsDNA Antibody ELISA 96 Determinations

INTENDED USE

An enzyme linked immunosorbent assay (ELISA) for the qualitative or semi-quantitative detection of IgG antibodies to double stranded DNA (dsDNA) in human serum, as an aid in the diagnosis of systemic lupus erythematosus (SLE) in conjunction with other laboratory and clinical findings.

SUMMARY AND EXPLANATION

Antinuclear antibodies (ANA) are found in various autoimmune diseases. ANA include antibodies to antigens of the nucleus such as to DNA, histone and various extractable nuclear antigens such as RNP, Sm, SS-A and SS-B. Three specificities occur with anti-DNA antibodies.

These include:

1. anti-dsDNA antibodies that react only with dsDNA
2. anti-ssDNA antibodies that react with ssDNA
3. anti-ds/ssDNA antibodies that react with both dsDNA and ssDNA.

Of these three types, anti-dsDNA antibodies are characteristic of systemic lupus erythematosus (SLE). They rarely occur in other autoimmune disorders.¹⁻⁶ The frequency and levels of these antibodies fluctuate with disease activity occurring overall in about 50-55% of SLE cases and in about 89% of SLE patients with active renal disease.³⁻⁷ Antibodies to dsDNA may disappear with immunosuppressive treatment and during remission. There is a good correlation between disease activity and anti-dsDNA antibody levels.⁸ The ImmuliTM dsDNA Antibody ELISA assay detects dsDNA antibodies of the IgG class. Anti-DNA antibodies of IgM and IgA isotypes also occur, but IgG class antibodies have been shown to be clinically relevant.³ Results are reported in International Units (IU)/ml. Both the calibrators and positive control have been calibrated against the World Health Organization (WHO) Reference Reagent Wo/80.⁹

PRINCIPLES OF PROCEDURE

The test is performed as a solid phase immunoassay. Microwells are coated with purified dsDNA antigen followed by a blocking step to reduce non-specific binding during the assay run. Controls, calibrators and patient sera are incubated in the antigen coated wells to allow specific antibodies present in the serum to bind to the dsDNA antigen. Unbound antibodies and other serum proteins are removed by washing the microwells. Bound antibodies are detected by adding an enzyme labeled anti-human IgG conjugate to the microwells. Unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the

EN

wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in International Units per milliliter (IU/ml) and reported as positive or negative.

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.¹⁰

Stop Solution is a dilute sulfuric acid solution. Sulfuric acid (H₂SO₄) is poisonous and corrosive. Do not ingest and avoid contact with skin and eyes. Avoid exposure to bases, metals or other compounds that may react with acids.

TMB Enzyme Substrate contains an irritant that may be harmful if inhaled, ingested or absorbed through the skin. Do not ingest and avoid contact with skin and eyes.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

Materials provided

ImmuliSTM dsDNA Antibody ELISA

REF 5120

Kits contains sufficient reagents to perform 96 determinations.

12 x 8

MICROPLATE dsDNA

Microplate with individual breakaway microwells. Coated with dsDNA antigen. Ready for use.

1 x 1.75 ml

CONTROL + dsDNA

Ready to use **Positive Control** (red cap). Contains human serum positive for dsDNA antibodies. The expected concentration range in IU/ml is printed on the label.

EN

1 x 1.75 ml		Ready to use Negative Control (<i>white cap</i>). Contains human serum.
5 x 1.75 ml	 	Ready to use set of 5 Calibrators. Calibrator A (green cap) 400 IU/ml, Calibrator B (violet cap) 200 IU/ml, Calibrator C (blue cap) 100 IU/ml, Calibrator D (yellow cap) 50 IU/ml, and Calibrator E (orange cap) 1 IU/ml. Derived from human serum containing dsDNA antibodies. Concentrations in IU/ml are printed on the labels.
1 x 15 ml		HRP goat anti-human IgG Conjugate. Ready for use. Color coded pink.
1 x 60 ml		Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml		TMB enzyme substrate. Ready for use. Protect from light.
1 x 15 ml		Stop Solution. Ready for use.
2 x vials		Powder Wash Buffer. Reconstitute to one liter each.
1 x		Protocol Sheets

Optional Components

1 x 60ml		Liquid concentrated Wash Buffer. Reconstitute to one liter.
----------	--	--

Symbols used on labels

	Lot number
	Catalog number
	In vitro diagnostic use
	Use by
	Storage temperature
	Read instructions before use
	Number of tests
	Manufacturer

Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels

EN

- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl per well.

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year.

PROCEDURE

Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- ***Good washing technique is critical.*** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. ***An automated microplate washer is recommended.***
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

Test Method

- Step 1** Let all reagents and specimens equilibrate to room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.
- or**

EN

For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

Qualitative				Semi-Quantitative			
A	Blank	S5	Etc.	A	Blank	S1	Etc.
B	-	S6		B	-	S2	
C	+	S7		C	+	S3	
D	Cal D	S8		D	Cal A	S4	
E	S1	S9		E	Cal B	S5	
F	S2	S10		F	Cal C	S6	
G	S3	S11		G	Cal D	S7	
H	S4	S12		H	Cal E	S8	
	1	2	3		1	2	3

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500ul** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.
- Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

EN

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <30 IU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining IU/ml. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

RESULTS

Calculations

The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

Abs. of Test Sample

----- X IU/ml of Calibrator D = IU/ml Test Sample

Abs. of Calibrator D

It is recommended that qualitative results be reported as “positive” or “negative.” Sample results greater than or equal to Calibrator D are considered positive.

2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in IU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as “positive,” “negative,” or “indeterminate” with IU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods, as indicated in **Limitations of Procedure**.

Interpretation

Interpretation values were determined by testing 181 normal blood donors and non-SLE disease control specimens. The mean of the normal subjects plus 3 SD was established as the assay cutoff and assigned a value of 50 IU/ml. IMMCO suggests use of the reference range below. Each laboratory should validate assay values for their own conditions.

<u>dsDNA Ab value</u>	<u>Interpretation</u>
<50 IU/ml	Negative
50-60 IU/ml	Indeterminate (Borderline)
>60 IU/ml	Positive

EN

Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining IU/ml values, multiply the units obtained by the dilution factor.

LIMITATIONS OF PROCEDURE

CLIA Complexity: High. CDC Analyte Identification Code: 0425. CDC Test System Identification Code: 28276.

The assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only. When the results of testing are indeterminate, further testing for the following is suggested:

- Antinuclear Antibodies on HEp-2 cells
- Antinuclear Antibodies on Mouse Kidney or Mouse Liver tissue sections
- Anti-ENA: RNP, Sm, SS-A (Ro), SS-B (La) Antibody ELISAs
- nDNA by indirect immunofluorescence

Above tests are available from IMMCO Diagnostics. Refer to Product Catalog.

Testing for complement levels C3 and C4, CH50 and immune complexes is also suggested. Strong positive results are indicative of SLE. However, negative results cannot necessarily rule out a diagnosis of SLE. When there is a high suspicion of SLE, then other tests, such as for ANA, ENA, nDNA and complement levels should also be considered. The results obtained serve only as an aid in the diagnosis. Taken alone, these results should not be interpreted as diagnostic.

EXPECTED VALUES AND SUMMARY OF CLINICAL STUDIES

The incidence of dsDNA antibodies varies depending upon the patient population. A compilation of incidence taken from the literature appears below.^{3,7,11-13}

Disease	% Incidence
Systemic Lupus Erythematosus	50-55
Active renal disease	89
Active non-renal disease	56
Inactive disease	30
Possible SLE	32
Rheumatoid Arthritis	<12
Systemic Scleroderma	<12
Normal patients	<1

Sets of clinical samples were tested on the ImmuLisa™ dsDNA Antibody ELISA. Results demonstrating incidence in the populations are provided below:

EN

Diagnosis	No. Tested	No. Positive	% Positive*
SLE	249	155	62.2
Primary anti-phospholipid syndrome	25	0	0
Primary myositis	20	0	0
Scleroderma	4	1	25.0
Rheumatoid arthritis	43	5	11.6
Graves' disease	9	0	0
Hashimoto's thyroiditis	14	1	7.1
Undifferentiated thyroid disease	9	1	11.1
Autoimmune hearing loss	8	0	0
Normal	159	6	3.8

* Indeterminate/borderline specimens considered positive

PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ dsDNA Antibody ELISA was evaluated by testing well-characterized specimens from SLE subjects with positive nDNA IFA results alongside disease controls and “normal” human sera. These specimens were also tested on commercially available ELISA kits. Only specimens in the linear range of the assay were included in the method comparison. These results are summarized below.

Immulisa™ dsDNA Antibody ELISA vs. other dsDNA antibody ELISA:

		Other dsDNA Antibody ELISA		
		Positive	Negative	Total
IMMCO	Positive*	71	10	81
dsDNA Ab.	Negative	6	86	92
ELISA	Total	77	96	173
Positive Percent Agreement:		92.2% (95% CI 83.2%-96.8%)		
Negative Percent Agreement:		89.6% (95% CI 81.3%-94.6%)		
Overall Percent Agreement:		90.8% (95% CI 85.2%-94.5%)		

* Indeterminate/borderline samples were considered positive.

SLE Subjects: 87

Disease Controls: 28

Healthy Normal Subjects: 58

Precision

Precision was tested with 6 specimens selected throughout the range of the assay. Runs of 6 replicates of each specimen were conducted in 12 assays. Repeatability was determined with 12 replicates of each specimen.

Sample	Mean (IU/ml)	Total Imprecision		Between days		Within run (Repeatability)	
		SD (IU/ml)	CV%	SD (IU/ml)	CV%	SD (IU/ml)	CV%
1	33.6	3.87	11.5%	4.04	11.7%	2.71	7.9%
2	41.0	4.43	10.8%	4.03	10.1%	4.08	10.2%
3	59.8	3.81	6.4%	3.83	6.6%	3.68	6.3%
4	124.9	10.31	8.3%	10.76	8.5%	5.56	4.4%
5	270.7	16.12	6.0%	16.44	5.9%	14.39	5.2%
6	361.6	24.99	6.9%	26.38	7.3%	14.72	4.1%

EN

Reproducibility

Assays of samples in the low negative range, moderate positive range and approximately +/- 20% of the assay cutoffs were performed to determine qualitative reproducibility. 86 replicates of each sample were tested. Assay results for the low negative and moderate positive specimens produced 100% qualitative agreement. Assay results for the cutoff -20% specimen produced 96.4% qualitative (negative) agreement. Assay results for the cutoff +20% specimen produced 98.8% qualitative (positive) agreement.

Linearity and Recovery

Studies were performed using equidistant dilution series of positive samples with values throughout the calibrator range to determine linear range of the assay. The linear range of the assay was determined to be from 13.6 – 450 IU/ml. Results are summarized below.

Sample	Test Range	Slope (95% CI)	Y-Intercept (95% CI)	R ²	% Recovery
1	13.7 - 115.0	0.980 (0.940 - 1.019)	0.171 (-2.699 - 3.042)	1.00	96.3 - 107.1
2	14.9 - 242.3	1.031 (0.983 - 1.079)	0.835 (-6.000 - 7.667)	1.00	91.2 - 106.3
3	43.4 - 398.3	0.911 (0.874 - 0.949)	4.088 (-5.331 - 13.507)	0.998	100.0 - 112.3
4	41.7 - 332.2	0.980 (0.915 - 1.045)	9.735 (-3.518 - 22.988)	0.993	87.6 - 110.6
5	48.9 - 399.5	0.911 (0.868 - 0.954)	6.865 (-4.132 - 17.862)	0.997	97.6 - 110.8

Limit of Detection

The limit of detection (LoD) for dsDNA antibody using this assay was determined to be 13.6 IU/ml based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples.

Interference

Interference was studied by mixing sera with known dsDNA antibody levels with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342 µmol/L), and Rheumatoid Factor (100 EU/ml).

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REFERENCES

1. Feltkamp TEW (Ed). The significance of the determination of anti-DNA and DNA/anti-DNA complexes. *Scand J Rheumatol Suppl.* 1975; 11:7-64.
2. Tan EM. Antinuclear antibodies: diagnostic markers for autoimmune diseases and probes for cell biology. *Adv Immunol.* 1989; 44:93-151.
3. Tzioufas AG, Manoussakis MN, Drosos AA et al. Enzyme immunoassays for the detection of IgG and IgM anti-dsDNA antibodies: clinical significance and specificity. *Clin Exp Rheumatol.* 1987; 5:247-253.
4. Ruffatti A, Calligaro A, Ross D et al. Anti-double stranded DNA antibodies in the healthy elderly: prevalence and characteristics. *J Clin Immunol.* 1990; 10:300-303.
5. Kadlubowski M, Jackson M, Yap PL and Neill G. Lack of specificity for antibodies to double stranded DNA found in four commercial kits. *J Clin Path.* 1991; 44:246-250.
6. Brinkman K, Termaat R, Van den Brink H et al. The specificity of the anti-dsDNA ELISA: a closer look. *J Immunol Methods.* 1991; 13:91-100.
7. Lange A. Evaluation of the simultaneous estimation of anti-dsDNA and anti-ssDNA antibodies for clinical purposes. *Clin Exp Immunol.* 1978; 31: 472-481.
8. Borg EJ, Horst G, Hum EJ et al. Measurement of increases in anti-double stranded DNA antibody level as a predictor of disease exacerbation in systemic lupus erythematosus: a long term, prospective study. *Arth Rheum.* 1990; 33:634-643.
9. Feltcamp TE, Kirkwood TB et al. The first international standard for antibodies to double stranded DNA. *Ann Rheum Dis.* 1988; 47:740-746.
10. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health. 2007; (HHS Pub No [CDC] 93-8395).
11. Conrad K, Schössler W, Hiepe F. Autoantibodies in systemic autoimmune diseases, a diagnostic reference. 2002; 69-72.
12. Vánca A et al. Myositis-specific and myositis-associated antibodies in overlap myositis in comparison to primary dermatomyositis: Relevance for clinical classification: retrospective study of 169 patients. *Joint Bone Spine.* 2010; 77:125-30.
13. Ruffatti A et al. Anti-double-stranded DNA antibodies in the healthy elderly: prevalence and characteristics. *J Clin Immunol.* 1990; 10:300-3.

For technical assistance please contact:

IMMCO Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immcodiagnostics.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante Autorizzato/Representante
Autorizado/Représentant Autorisé
EMERGO Group, Inc.
Molenstraat 15, 2513 BH, The Hague,
The Netherlands
Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com

OCT2011

Document No. PI5120